

Product Name (必要項目)	Description (必要項目)	Caution (必要項目)	Catalog No. (必要項目)	Model (必要項目)
Sanplatec Dry Keeper	<ul style="list-style-type: none"> This is one of the most popular models. The electronic dehumidification system. The dehumidification system. 1 set comes with 3 adjustable shelves. Large, in-door digital hygrometer. The shelf can be repositioned. 	<ul style="list-style-type: none"> Solid polymer electrolyte membrane. Do not store volatile chemicals. 	00001	Auto C-3B
Sanplatec Dry Keeper			00002	Auto C-3W
Sanplatec Dry Keeper			00003	Auto C-3BS
Sanplatec Dry Keeper			00004	Shelf only(C-3)
Sanplatec Dry Keeper	<ul style="list-style-type: none"> This is one of the most popular models. The electronic dehumidification system. The dehumidification system. 1 set comes with 3 adjustable shelves. Additional shelves and rails are available. Large, in-door digital hygrometer. The shelf can be repositioned. 	<ul style="list-style-type: none"> Solid polymer electrolyte membrane. Do not store volatile chemicals. 	00010	Auto A-3B
Sanplatec Dry Keeper			00011	Auto A-3W
Sanplatec Dry Keeper			00012	Shelf (A-3)
Sanplatec Dry Keeper	<ul style="list-style-type: none"> This is one of the most popular models. 1 set comes with 3 adjustable shelves. Large, in-door digital hygrometer. The shelf can be repositioned. 		00030	Non-auto C-3B
Sanplatec Dry Keeper			00031	Non-auto C-3W
Sanplatec Dry Keeper			00032	Non-auto C-3BS
Sanplatec Dry Keeper			00004	Shelf only(C-3)
Sanplatec Dry Keeper	<ul style="list-style-type: none"> This is one of the most popular models. 1 set comes with 3 adjustable shelves. Large, in-door digital hygrometer. The shelf can be repositioned. 		00040	A-3B
Sanplatec Dry Keeper			00041	A-3W
Sanplatec Mini Desiccator (Mini Desi)	<ul style="list-style-type: none"> Compact and light. Example application: Carry over. Each type has a silica gel. We can make 2 or 3 steps. 		00050	MD-1
Sanplatec Mini Desiccator (Mini Desi)			00051	MD-1H
Sanplatec Mini Desiccator (Mini Desi)	<ul style="list-style-type: none"> Smoke and transparent type. Performance is the same as smoke type. Compact and light. Example application: Carry over. Each type has a silica gel. We can make 2 or 3 steps. 		00056	Smoke type MD-1S
Sanplatec Mini Desiccator (Mini Desi)			00057	Smoke type MD-1SH
Polycarbonate Vacuum Desiccator			95001	PC-150・210・250
Polycarbonate Vacuum Desiccator	<ul style="list-style-type: none"> Clear view. Easy to see the contents. Dry seal type with O-ring. Wet seal type. The highest intensity in transparent. The ceramic plate has heat resistance. PC material feature: good purification. Spare plate for each type (PC-210K). Nozzle diameter is for 9mm. 		00062	PC-210K
Polycarbonate Vacuum Desiccator			00063	PC-250K
Polycarbonate Vacuum Desiccator			00064	PC-210
Polycarbonate Vacuum Desiccator			00065	PC-250
Polycarbonate Vacuum Desiccator			00066	Plate only 210 type
Polycarbonate Vacuum Desiccator			00067	Plate only 250 type
Polycarbonate Vacuum Desiccator			04705	PC-210KG
Polycarbonate Vacuum Desiccator			04706	PC-250KG
Middle Size Desiccator	<ul style="list-style-type: none"> Adjustable shelf at 28mm interval. The door has a key. It is safe. 		00118	Standard C-4
Middle Size Desiccator			00119	Standard C-5
Middle Size Desiccator			00120	Standard C-6
Middle Size Desiccator			00121	Shelf (C-4 and 5 type)
Middle Size Desiccator			00122	Shelf (C-6 type)

Middle Size Auto Desiccator	<ul style="list-style-type: none"> ●Full-automation and no need ●It maintains RH levels at 30 		00128	Auto C-4
Middle Size Auto Desiccator			00129	Auto C-5
Middle Size Auto Desiccator			00130	Auto C-6
Middle Size Non-Auto Desiccator	<ul style="list-style-type: none"> ●The shelf can be repositione ●High-cost performance & hig 		00131	AC-500type
Middle Size Non-Auto Desiccator			00132	AC-1000type
D-BOX (Auto type)	<ul style="list-style-type: none"> ●Dehumidification system with ●Compact and slim body ●open and close the door with ●Stackable ●The main body is manufactu 		00021	Shelf (D-BOX·i-BOX)
D-BOX (Auto type)			00150	ver.1
i-BOX (Auto type)	<ul style="list-style-type: none"> ●Dehumidification system with ●It maintains relative humidity ●Stackable 	●Please select your preferred c	00151	Forest
i-BOX (Auto type)			00152	Tropical
i-BOX (Auto type)			00153	Ocean
i-BOX (Auto type)			00154	Snow
i-BOX (Auto type)			00155	Moon
P-BOX	<ul style="list-style-type: none"> ●Made in Japan with high-qua ●Compact size even when pla ●Lightweight 5kg for women-f ●The back of the main body h ●Application: culture, reagent 		00156	Vertical type
P-BOX			00157	Side type
P-BOX			00158	Shelf only (Vertical type)
P-BOX			00159	Shelf only (Side type)
Acrylic Vacuum Desiccator	<ul style="list-style-type: none"> ●This is a vacuum desiccator ●It is also usable as a package ●Vacuum gauge 1pcs for 0~ 		00175	UB type
Acrylic Vacuum Desiccator	<ul style="list-style-type: none"> ●This is a vacuum desiccator ●It is also usable as a package ●Vacuum gauge 1pcs for 0~ 		00176	MB type
Acrylic Vacuum Desiccator	<ul style="list-style-type: none"> ●This is a vacuum desiccator ●Vacuum gauge 1 pcs for 0~ 		00177	YB-1 type
Acrylic Vacuum Desiccator			00178	YB-2 type
Acrylic Vacuum Desiccator			00179	YB-3 type
Acrylic Vacuum Desiccator			00180	shelf for YB type
Acrylic Vacuum Desiccator	<ul style="list-style-type: none"> ●This is a vacuum desiccator ●Vacuum gauge 1pcs for 0~ 		00181	SB type
Nonporous Package Tester	<ul style="list-style-type: none"> ●This product provides a simp ●The transparent test chamb ●Precise pressure adjustment ●The lamp display allows you ●Connect to your vacuum pur ●Chambers can be selected a ●The timer allows you to set ●Spec:AC100V·50/60Hz ●Set details (Code:00182): co 		00182	Control Box
Nonporous Package Tester			00183	Chamber type 1
Nonporous Package Tester			00184	Chamber type 2
Nonporous Package Tester			00185	Chamber type 3
Glove Box	<ul style="list-style-type: none"> ●The opening part is large to put in and out samples or equipment. ●With a lock ●Accessory: 1pair of gloves 		00186	ST-1 type
Glove Box			00187	ST-2 type
Glove (Latex)	<ul style="list-style-type: none"> ●Latex gloves 		30226	7inch
Glove (Latex)			30227	8inch
Glove Box Accessory	<ul style="list-style-type: none"> ●Made from latex 		00193	6 inch
Glove Box Accessory			00194	ST-1 Packing
Glove Box Accessory			00195	ST-2·DX Packing

Glove Box Accessory			00196	Packing for side box
Acrylic Vacuum Glove Box, Clear View	<ul style="list-style-type: none"> ● This glove box allows two people to work at the same time. ● Work can be performed in an inert gas atmosphere. ● Work equipment can be moved in and out from the left door. Instruments and samples can be moved in and out of main body without changing the environment by using pass box on the right side of the door. ● Both main body and pass box are designed to be used safely for vacuum-pumping. ● Accessories: Glove 2pairs, 4: O-rings for glove retaining, 2 valves for air supply, 2 valves for exhaust, 2 valves for glove port, and vacuum gauge 2pcs. 	<ul style="list-style-type: none"> ● Do not operate when the inside is under decompression. 	00199	400 type
Mini Hood	<ul style="list-style-type: none"> ● Compact and convenient one. ● It is made of vinyl chloride acrylate. ● The front door is slidable and lockable. ● A power outlet (Japanese Standard). ● Spec: Fan:AC100V 50/60Hz 	<ul style="list-style-type: none"> ● As metal is used partially such as door, Please replace parts when a metal part is damaged. ● Do not use flammable materials. 	00202	MF-S type
Mini Hood			00203	MF-L type
Acrylic Column	<ul style="list-style-type: none"> ● Parts detail: main body 1pc (No.1), Perforated plate 2pcs, Gasket 1pc. 		00205	No.1
Acrylic Column			00206	No.2
Acrylic Column			00207	No.3
Acrylic Column			00208	No.4
Acrylic Column			00209	No.5
Acrylic Column			00210	No.6
Acrylic Column			00211	No.7
Acrylic Column			00212	No.8
Acrylic Vacuum Desiccator	<ul style="list-style-type: none"> ● This is a vacuum desiccator made of acrylic made by hand with advanced bonding techniques. ● It is also usable as a package tester because it is designed to depressurize even when liquid is inside. ● Vacuum gauge 1pc for 0~ -0.1Mpa, Valve 2pcs, nozzle O.D. 9mm, Gasket: Silicone ● Suitable to store the moisture sensitive materials. ● Perfect compact design for storage and management in various parts such as sampling, sensors, optical jigs, and gauges. Smaller types are also available. Please contact us. ● Applications: Available in a wide range, including vacuum storage, testing, and manage gas purge sampling. 		00225	SNS type
Acrylic Vacuum Desiccator			00226	SNH type
Acrylic Vacuum Desiccator			00227	SNW type
Acrylic Vacuum Desiccator			04710	SNH-TAP type
Acrylic Vacuum Desiccator			04711	SNW-TAP type
Glove Box	<ul style="list-style-type: none"> ● Sterilization lamp plates are available. ● Usable as aseptic box and decontamination box. ● Large opening ensures safe work. ● Accessory: Glove 1pair, Air e ● Germicidal lamp, Fluorescent lamp ● Pass box: (Inner dimension) ● Fluorescent lamp: 10W ● UV: 10W ● Fan: 50/60Hz ● Please contact us for detail. 		00228	U-DX type

Acrylic Vacuum Desiccator (Portable type)	<ul style="list-style-type: none"> ● This is an acrylic vacuum desiccator ● The handles for easy carrying ● Vacuum gauge 1pc (for 0~ -0.1mmHg) 		00230	SB-II P type
Polycarbonate Desiccator Mini	<ul style="list-style-type: none"> ● With Spigot type : nozzle diameter 10mm ● With vacuum gauge type included ● Each type has a thick performance ● Dry seal type with O-ring. 		00240	PC-150K
Polycarbonate Desiccator Mini			00241	PC-150
Polycarbonate Desiccator Mini			04707	PC-150KG
Table Top Clean Bench (safety sliding door)	<ul style="list-style-type: none"> ● The sliding door is adjustable at 28mm intervals. ● The clean unit blows clean air, removing more than 99.99% of the dust and debris on 0.3 μm suspended in the atmosphere. ● Frame: aluminum ● PVC clear: door, side, top, back ● PVC white: bottom board 		00242	500 type
Table Top Clean Bench (safety sliding door)			00243	700 type
Table Top Clean Bench (safety sliding door)			00244	1000 type
Table Top Clean Bench (safety sliding door)			00245	1000H type
Table Top Hood	<ul style="list-style-type: none"> ● The sliding door is adjustable ● Can be installed on a laboratory ● Spec: aluminum frame, PVC ● PVC white: bottom board ● AC fan: 2.8m³ /min ● Accessory: 100 φ × 3m duct 	● As metal is used partially such as	00246	500H type
Table Top Hood			00247	700 type
Table Top Hood			00248	1000 type
Table Top Hood			00249	1000H type
Compact Desiccator	<ul style="list-style-type: none"> ● One-touch door type for easy operation ● Adjustable shelves ● Water-washable ● Accessory: silica gel and 2 pieces ● The color of main body: gray 		00250	HA4-B01
Compact Desiccator			00251	HA5-B81
Acrylic Vacuum Desiccator (pipe type)	<ul style="list-style-type: none"> ● Perfect vacuum desiccator type ● Vacuum gauge 1pc (for 0~ -0.1mmHg) 		00252	PB type
Acrylic Vacuum Desiccator (pipe type)			00392	PS type
Acrylic Vacuum Desiccator (pipe type)			04712	PS-II type
Clean Bench (antistatic sheet)	<ul style="list-style-type: none"> ● Simplified Clean Bench with ● Easy setting ● Filter: HEPA filter, Non-woven ● Airflow: 1m³ /min ● Power :AC100V ● Frame: aluminum(antistatic type) ● Cover: 0.3t antistatic vinyl sheet 		00263	S-800
Clean Bench (antistatic sheet)			00264	S-1200
Clean Bench (PVC antistatic board)	<ul style="list-style-type: none"> ● Clean Bench incorporating PVC ● Filter: HEPA filter, Non-woven ● Airflow: 1m³ /min ● Power: AC100V ● Frame: aluminum(antistatic type) ● Plate: PVC Antistatic 3mm(s) 		00265	P-800
Clean Bench (PVC antistatic board)			00266	P-1200
E-BOX (Non-auto type)	<ul style="list-style-type: none"> ● Easy storage and compact. 	● E-BOX does not include an a	00315	Standard
i-BOX (Non-auto type)	<ul style="list-style-type: none"> ● Non-Auto type desiccator. Forest ● Flat top surface. Stackable u ● Easy to open and close with ● Please choose your favorite 		00316	Forest
i-BOX (Non-auto type)			00317	Tropical
i-BOX (Non-auto type)			00318	Ocean
i-BOX (Non-auto type)			00319	Snow
i-BOX (Non-auto type)			00320	Moon
Acrylic Column		● Be cautious not to put excess	00340	No. 1
Acrylic Column			00341	No.2

Acrylic Column			00342	No.3
Acrylic Column			00343	No.4
Acrylic Column			00344	No.5
Acrylic Column			00345	No.6
Large Size Desiccator ALW type	<ul style="list-style-type: none"> ● This desiccator allows quick ● The shelf is made of stainless ● The aluminum outer frame is ● Auto type is equipped with d ● Maintain the inside humidity ● Door: Magnetic type 	<ul style="list-style-type: none"> ● Gas replacable type does not 	00273	ALW-G-NK
Large Size Desiccator ALW type			00351	ALW-X-NK
Large Size Desiccator ALW type			00371	ALW-A-NK
Large Size Desiccator ALS type	<ul style="list-style-type: none"> ● This desiccator allows quick ● The shelf is made of stainless ● The aluminum outer frame is ● Auto type is equipped with d ● Maintain the inside humidity ● Door: Magnetic type 	<ul style="list-style-type: none"> ● Gas replacable type does not 	00274	ALS-G-NK
Large Size Desiccator ALS type			00352	ALS-X-NK
Large Size Desiccator ALS type			00370	Shelf only (AL type)
Large Size Desiccator ALS type			00372	ALS-A-NK
Middle Size Desiccator ALG-1 type	<ul style="list-style-type: none"> ● This desiccator allows quick ● The shelf is made of stainless ● The aluminum outer frame is ● Auto type is equipped with d ● Maintain the inside humidity ● Door: Magnetic type 		00275	ALG1-G-N
Middle Size Desiccator ALG-1 type			00357	ALG1-X-N
Middle Size Desiccator ALG-1 type			00358	ALG1-X-NK
Middle Size Desiccator ALG-1 type			00373	ALG1-A-N
Middle Size Desiccator ALG-3 type	<ul style="list-style-type: none"> ● This desiccator allows quick ● The shelf is made of stainless ● The aluminum outer frame is ● Auto type is equipped with d ● Maintain the inside humidity ● Door: Magnetic type 		00277	ALG3-G-N
Middle Size Desiccator ALG-3 type			00360	ALG3-X-N
Middle Size Desiccator ALG-3 type			00361	ALG3-X-NK
Middle Size Desiccator ALG-3 type			00375	ALG3-A-N
Acrylic Vacuum Desiccator	<ul style="list-style-type: none"> ● Perfect vacuum desiccator to store moisture-sensitive materials. 		00390	SNO type
Acrylic Vacuum Desiccator			04713	SNO-TAP type
Acrylic Vacuum Desiccator	<ul style="list-style-type: none"> ● This is a vacuum desiccator ● Vacuum gauge 1pcs for 0~ 		00391	SB-II type
Acrylic Vacuum Desiccator			04714	SB-II-TAP type
Acrylic Vacuum Defoaming Device	<ul style="list-style-type: none"> ● The status of stirring and de 	<ul style="list-style-type: none"> ● Do not include a stirring moto ● The outer dimension is only th 	00395	SNO-VS type
Acrylic Vacuum Defoaming Device			00396	SNH-VS type
Acrylic Vacuum Defoaming Device			00397	SNW-VS type
Chemical Graduated Cylinder			01001	10mL
Chemical Graduated Cylinder			01002	20mL
Chemical Graduated Cylinder			01003	50mL
Chemical Graduated Cylinder			01004	100mL
Chemical Graduated Cylinder			01005	200mL
Chemical Graduated Cylinder			01006	300mL
Chemical Graduated Cylinder			01007	500mL
Chemical Graduated Cylinder			01008	1L
Chemical Graduated Cylinder			01009	2L

TPX® Graduated Cylinder	<ul style="list-style-type: none"> ●TPX graduated cylinder series ●Great chemical resistance, high strength ●The upper part of 10mL and 20mL ●The graduation can be also used ●The bottom stand part is for 	01011	10mL
TPX® Graduated Cylinder		01012	20mL
TPX® Graduated Cylinder		01013	50mL
TPX® Graduated Cylinder		01014	100mL
TPX® Graduated Cylinder		01015	200mL
TPX® Graduated Cylinder		01016	300mL
TPX® Graduated Cylinder		01017	500mL
TPX® Graduated Cylinder		01018	1L
TPX® Graduated Cylinder		01019	2L
PP Graduated Beaker	<ul style="list-style-type: none"> ●Fine graduations enable easy viewing. ●Graduate scale on a max scale is surrounded a cup. ●It can also be used when stirring with a magnetic stirrer. 	01021	100mL
PP Graduated Beaker		01022	200mL
PP Graduated Beaker		01023	300mL
PP Graduated Beaker		01024	500mL
PP Graduated Beaker		01025	1L
PP Graduated Beaker with Handle	<ul style="list-style-type: none"> ●The handle type is easy to use ●Fine graduations allow users ●Usable for agitation with magnetic stirrer 	01026	100mL
PP Graduated Beaker with Handle		01027	200mL
PP Graduated Beaker with Handle		01028	300mL
PP Graduated Beaker with Handle		01029	500mL
PP Graduated Beaker with Handle		01030	1L
PP Graduated Beaker with Handle		01031	2L
PP Graduated Beaker with Handle		01032	3L
PP Graduated Beaker with Handle		01033	5L
TPX® Graduated Beaker	<ul style="list-style-type: none"> ●Plastic integrally molded beaker ●Lightweight to handle easily ●Usable for agitation with a magnetic stirrer ●The heat-resistant temperature ●Autoclavable 	01041	100mL
TPX® Graduated Beaker		01042	200mL
TPX® Graduated Beaker		01043	300mL
TPX® Graduated Beaker		01044	500mL
TPX® Graduated Beaker		01045	1L
Manifold Filtration Device	<ul style="list-style-type: none"> ●Multiple samples can be filtered simultaneously. ●Two nozzles are attached for depressurization and drainage. ●The opposite side of the nozzle is a screw cap, which allows cleaning of the inside of the main body tube. 	01049	3 pipes
Manifold Filtration Device		01050	5 pipes
TPX® Graduated Beaker with Handle	<ul style="list-style-type: none"> ●Plastic integrally molded beaker ●Lightweight to handle easily ●Usable for agitation with a magnetic stirrer ●The heat-resistant temperature ●Autoclavable 	01051	100mL
TPX® Graduated Beaker with Handle		01052	200mL
TPX® Graduated Beaker with Handle		01053	300mL
TPX® Graduated Beaker with Handle		01054	500mL
TPX® Graduated Beaker with Handle		01055	1L

TPX® Graduated Beaker with Handle		01056	2L
TPX® Graduated Beaker with Handle		01057	3L
TPX® Graduated Beaker with Handle		01058	5L
Beaker Lid	●Lid for TPX® beaker, PP bea	01061	100mL
Beaker Lid		01062	200mL
Beaker Lid		01063	300mL
Beaker Lid		01064	500mL
Beaker Lid		01065	1L
Beaker Lid		01066	2L
96 Well Plate Cover	●Usable for multi-hole plates	01073	Non-sterilized (50pcs)
96 Well Plate Cover		01074	Sterilized (50pcs)
Square Dish	●Gamma sterilized ●Good stability when stacking ●Example applications: store various cell kinds (ex, agar cultures, floating cell cultures, microbial cultures, cells, microorganisms and reagents, etc.)	01075	MPS01F01S (1pc × 50)
Square Dish		01076	MPS01F04S (5pcs × 16)
Sanplatec Large Acrylic Graduated Cylinder	●The handle is convenient for	01077	3L
Sanplatec Large Acrylic Graduated Cylinder		01078	5L
Sanplatec Large Acrylic Graduated Cylinder		01079	10L
Graduated Cylinder with Spigot		01080	1L
Graduated Cylinder with Spigot		01081	2L
Graduated Cylinder with Spigot		01082	3L
Graduated Cylinder with Spigot		01083	5L
Graduated Cylinder with Spigot		01084	10L
TPX® Conical Centrifuge Tube	●It can be centrifuged up to 2 ●The tip of the tube is gradua ●Resistance to organic solven	01089	15mL(300pcs)
TPX® High Speed Funnel	●The heat resistance tempera ●The inside is formed with he ●It is excellent in chemical re ●Autoclavable	01091	60mm
TPX® High Speed Funnel		01092	75mm
TPX® High Speed Funnel		01093	90mm
TPX® High Speed Funnel		01094	120mm
TPX® High Speed Funnel		01095	150mm
TPX® High Speed Funnel		01096	180mm
TPX® High Speed Funnel		01097	210mm
TPX® High Speed Funnel		01098	240mm
TPX® High Speed Funnel		01099	300mm
PP High Speed Funnel	●The heat resistance temperature: 120°C ●The inside is formed with helical ribs to prevent the filter from sticking.	01101	60mm
PP High Speed Funnel		01102	75mm
PP High Speed Funnel		01103	90mm
PP High Speed Funnel		01104	120mm
PP High Speed Funnel		01105	150mm
PP High Speed Funnel		01106	180mm
PP High Speed Funnel		01107	210mm
PP High Speed Funnel		01108	240mm

PP High Speed Funnel			01109	300mm
Funnel Stand	<ul style="list-style-type: none"> ●The size of the blade part is ●The base part is made of PVC 		01126	KVR-1 type
Funnel Stand			01127	KVR-2 type
Funnel Stand			01128	KVR-5 type
Funnel Stand	<ul style="list-style-type: none"> ●The size of the blade part is common for each type and exchangeable. ●The base part is made of PVC. 		01131	Y-1 type
Funnel Stand			01132	Y-2 type
Funnel Stand			01133	Y-3 type
Funnel Stand			01134	Y-4 type
Funnel Stand			01135	Y-5 type
Funnel Stand	<ul style="list-style-type: none"> ●2pcs (diameter :36mm) of ite 		01136	Standard
Drying Rack for Lab Equipment	<ul style="list-style-type: none"> ●Dry rack for equipment. Wall hanging type. ●Adjustable hangers for glass equipment size or qty. ●Hanger is sturdy (made of ABS resin) ●Add preferred qty of hangers if necessary. 		01155	S type
Drying Rack for Lab Equipment			01157	Hanger only
Drying Rack for Lab Equipment	<ul style="list-style-type: none"> ●Easy lateral motion of sink o ●It can be used for both verti 		01160	Assembly type
Graduated Pipette			01165	1mL
Graduated Pipette			01166	2mL
Graduated Pipette			01167	5mL
Graduated Pipette			01168	10mL
Graduated Pipette			01169	20mL
Volumetric Pipette			01170	5mL
Volumetric Pipette			01171	10mL
Volumetric Pipette			01172	20mL
Volumetric Pipette			01173	25mL
Volumetric Pipette			01174	50mL
Volumetric Pipette			01175	100mL
Pipette Stand	<ul style="list-style-type: none"> ●The height of blade is adjust ●The distance between 2 blade 		01176	S type
Bellows Pipette			01180	1mL(10pcs)
Bellows Pipette			01181	2mL(10pcs)
Bellows Pipette			01182	5mL(10pcs)
Bellows Pipette			01183	10mL(10pcs)
Bellows Pipette			01184	20mL(5pcs)
Bellows Pipette			01185	50mL(5pcs)
Volumetric Pipette		<ul style="list-style-type: none"> ●Silicon cap is necessary when 	01190	1mL
Volumetric Pipette			01191	2mL
Volumetric Pipette			01192	5mL
Volumetric Pipette			01193	10mL
Volumetric Pipette			01194	20mL
Pipette Washer (main body)	<ul style="list-style-type: none"> ●The Siphon part includes joint ●Transparent allows clear-vie 		01206	A-1
Pipette Washer (main body)			01207	A-2
Pipette Washer (main body)			01208	A-3
Pipette Washer (basket)		<ul style="list-style-type: none"> ●The dimension of the basket is the external dimension of the handle. 	01210	B-1
Pipette Washer (basket)			01211	B-2

Pipette Washer (basket)			01212	B-3
Pipette Washer (basket)			01213	B-4
Pipette Washer (basket)			01214	B-5
Pipette Washer (basket)	<ul style="list-style-type: none"> ●The Siphon part includes a jet ●Transparent enables clear-v 		01215	AB-1
Pipette Washer (basket)			01216	AB-2
Pipette Washer (basket)			01217	AB-3
PVC Tank	<ul style="list-style-type: none"> ●Usable for washing lab equip 	<ul style="list-style-type: none"> ●D-1/D-2: Gray 	01220	D-1
PVC Tank			01221	D-2
PVC Tank			01222	D-3
PVC Tank			01223	D-4
PVC Tank			01224	D-5
Tongs (PVC)	<ul style="list-style-type: none"> ●Vinyl chloride tong, acid resi 		01237	PVC type
Burette	<ul style="list-style-type: none"> ●The stopcock is made of PM 		01197	Standard
Burette			01240	10mL
Burette			01241	25mL
Burette			01242	50mL
Square PVC Washing Tank	<ul style="list-style-type: none"> ●Main body includes a lid. ●Basket with stirring rod ●All part is made of PVC. 		01246	SK-1 (Main Body)
Square PVC Washing Tank			01247	SK-2(Main Body)
Square PVC Washing Tank			01248	SB-1(Basket for SK-1)
Square PVC Washing Tank			01249	SB-2(Basket for SK-2)
Sanplatec PE Cup			01250	100mL (1000pcs)
Sanplatec PE Cup			01251	200mL (1000pcs)
Sanplatec PE Cup			01252	300mL (500pcs)
Sanplatec PE Cup			01253	500mL (200pcs)
Sanplatec PE Cup			01254	1L (200pcs)
PE Stirring Rod			01257	10MM
PE Stirring Rod			01258	12MM
Colorimetric Tube Stand	<ul style="list-style-type: none"> ●It comes with a black and w 		01260	50mL
Colorimetric Tube Stand			01261	100mL
Transparency Meter	<ul style="list-style-type: none"> ●Accessories: Drainage tube (150 mm), clamp (for drainage tube) ●Spec: transparency meter (cylinder (300mm, 500mm, 1000mm), frame, sign, tubing for drainage) ●Lighter and easier to handle than glass transparency meter 	<ul style="list-style-type: none"> ●Clamp color is either blue or green. 	01262	300 type (Set)
Transparency Meter			01263	500 type (Set)
Transparency Meter			01264	300 type (Cylinder)
Transparency Meter			01265	500 type (Cylindert)
Transparency Meter			01267	1000 type (Set)
Transparency Meter			01268	1000 type (Cylinder)
Mouse Restrainer with Defecation and Urination Port	<ul style="list-style-type: none"> ●Defecation / urination port is placed on the lower plate near the tail of the conventional adjuster ●Users are allowed to use ●Feces and urine prolapse as defensive response can be removed outside of equipment allowing to keep inside clean to use continuously. 	<ul style="list-style-type: none"> ●When cleaning, wipe out using a cloth impregnated with neutral detergent and rinse with water. ●Be aware that cleaning and disinfection with alcohol such as ethanol and other organic solvents may cause material deterioration. 	01318	Mause HKA

Sanpla Cup (PP Beaker)	<ul style="list-style-type: none"> ● Made of PP. Excellent chemical resistance ● Convenient with pour-spout ● Graduation scales on each side ● Reusable and Stackable to save space 		01660c	100mL (500pcs)
Sanpla Cup (PP Beaker)			01661c	200mL (500pcs)
Sanpla Cup (PP Beaker)			01662c	300mL (250pcs)
Sanpla Cup (PP Beaker)			01663c	500mL (200pcs)
Sanpla Cup (PP Beaker)			01664c	1L (100pcs)
PP Disposable Cup (PP Beaker)	<ul style="list-style-type: none"> ● PP material. Great chemical and heat resistance ● Convenient for agitation and preparation of samples in laboratories and for the preparation and agitation of dyestuffs and paints ● Good transparency ● The body is thick and balanced for easy-grip. ● The 1 L and 2L sides are embossed to make them easier to hold. 		01665c	100mL (1000pcs)
PP Disposable Cup (PP Beaker)			01666c	200mL (1000pcs)
PP Disposable Cup (PP Beaker)			01667c	300mL (500pcs)
PP Disposable Cup (PP Beaker)			01668c	500mL (250pcs)
PP Disposable Cup (PP Beaker)			01669c	1L (100pcs)
PP Disposable Cup (PP Beaker)			01840c	2L (100pcs)
PP Disposable Cup (PP Beaker)			01841c	50mL (1000pcs)
Boxcel	<ul style="list-style-type: none"> ● Suitable for packaging tip-pec ● 2 stretchable films pack con ● Great transparency for better ● Protect contents from impact ● Applications: transportation ● Precision parts (watch and e ● Pharmaceuticals and delicat ● Laboratory devices and che 	● The content size is for refer	01670	NO.23
Boxcel			01671	NO.32
Boxcel			01672	NO.33
Boxcel			01674	NO.98
Boxcel			01676	NO.148
Boxcel			01677	NO.149
Boxcel			01678	NO.225
Boxcel			01680	NO.438
Boxcel			01681	NO.467
Boxcel			01682	NO.554
Boxcel			01683	NO.651
Boxcel			01684	No.1180
Boxcel			01685	NO.12127
Eye Drop Container (Eito)	<ul style="list-style-type: none"> ● Body, cap: white color 		01743	0 type, No cover (100pcs)
Eye Drop Container (Eito)			01744	1 type, No cover (100pcs)
Eye Drop Container (Eito)			01745	2 type, No cover (100pcs)
Eye Drop Container (Eito)			01746	3 type, No cover (100pcs)
Small Beaker (Yakuhai)			01788	10mL(PS)
Small Beaker (Yakuhai)			01789	20mL(PS) 2000pcs
Small Beaker (Yakuhai)			01790	20mL(PP) 2000pcs
Small Beaker (Yakuhai)			01791	50mL(PP) 2000pcs

Square Bottle		●Cap, inner seal: not included	01924	KH-100 (730pcs)
Square Bottle			01925	KH-150 (540pcs)
Square Bottle			01926	KH-200 (430pcs)
Square Bottle			01927	KH-300 (284pcs)
Square Bottle			01928	KH-500 (170pcs)
Dropping Pipette			01975	1mL(1000pcs)
Dropping Pipette			01977	3mL(500pcs)
Dropping Pipette			01978	10mL(500pcs)
Dropping Pipette		● Graduation scales in each 0.	01976	2mL(500pcs)
PE Disposable Cup (PE Beaker)		● Not sterilized	01990	100mL (1000pcs)
PE Disposable Cup (PE Beaker)		● Disposable cup made of PE	01991	200mL (1000pcs)
PE Disposable Cup (PE Beaker)		● Soft and easy to pour	01992	300mL (500pcs)
PE Disposable Cup (PE Beaker)			01993	500mL (250pcs)
PE Disposable Cup (PE Beaker)			01994	1L (100pcs)
PE Disposable Cup (PE Beaker)			01995	2L (100pcs)
Sanplatec Clear Wide Mouth Bottle	●No inner seal, Worry-free of ●Plastic plasticizer (DOA·DO ●Heat resistant temperature	●Each size is for reference only ●The height is measured with v	02011c	30mL (400pcs)
Sanplatec Clear Wide Mouth Bottle			02012c	50mL (400pcs)
Sanplatec Clear Wide Mouth Bottle			02013c	100mL (200pcs)
Sanplatec Clear Wide Mouth Bottle			02014c	250mL (100pcs)
Sanplatec Clear Wide Mouth Bottle			02015c	500mL (100pcs)
Sanplatec Clear Wide Mouth Bottle			02016c	1L (50pcs)
Sanplatec Clear Wide Mouth Bottle			02017c	2L (30pcs)
Sanplatec PP Narrow Mouth Bottle	●No inner seal ●A highly transparent grade of ●Great heat and chemical res ●For a wide range of industrie	●Each size is for reference only	02030c	30mL (400pcs)
Sanplatec PP Narrow Mouth Bottle			02031c	50mL (400pcs)
Sanplatec PP Narrow Mouth Bottle			02032c	100mL (200pcs)
Sanplatec PP Narrow Mouth Bottle			02033c	250mL (100pcs)
Sanplatec PP Narrow Mouth Bottle			02034c	500mL (100pcs)
Sanplatec PP Narrow Mouth Bottle			02035c	1L (50pcs)
Sanplatec PP Narrow Mouth Bottle			02036c	2L (30pcs)
Sanplatec PP Wide Mouth Bottle	●No inner seal ●A highly transparent grade of ●Great heat and chemical res ●For wide range of industries	●Each size is for reference only	02040c	20mL (400pcs)
Sanplatec PP Wide Mouth Bottle			02041c	30mL (400pcs)
Sanplatec PP Wide Mouth Bottle			02042c	50mL (400pcs)
Sanplatec PP Wide Mouth Bottle			02043c	100mL (200pcs)
Sanplatec PP Wide Mouth Bottle			02044c	250mL (100pcs)
Sanplatec PP Wide Mouth Bottle			02045c	500mL (100pcs)
Sanplatec PP Wide Mouth Bottle			02046c	1L (50pcs)
Sanplatec PP Wide Mouth Bottle			02047c	2L (30pcs)
Sanplatec PE Narrow Mouth Bottle	●No inner seal for 30mL~50 ●5L, 10L have a hole for a stri	●Each size is for reference only ●The height is measured with v	02060c	30mL (400pcs)
Sanplatec PE Narrow Mouth Bottle			02061c	50mL (400pcs)
Sanplatec PE Narrow Mouth Bottle			02062c	100mL (200pcs)
Sanplatec PE Narrow Mouth Bottle			02063c	250mL (100pcs)

Sanplatec PE Narrow Mouth Bottle			02064c	500mL (100pcs)
Sanplatec PE Narrow Mouth Bottle			02065c	1L (50pcs)
Sanplatec PE Narrow Mouth Bottle			02066c	2L (30pcs)
Sanplatec PE Narrow Mouth Bottle			02067c	3L (30pcs)
Sanplatec PE Narrow Mouth Bottle			02068c	5L (15pcs)
Sanplatec PE Narrow Mouth Bottle			02069c	10L (8pcs)
Sanplatec PE Wide Mouth Bottle	<ul style="list-style-type: none"> ●No inner seal for 20mL~50mL ●5L and 10L have a hole for a 	<ul style="list-style-type: none"> ●Each size is for reference only ●The height is measured with v 	02080c	20mL (400pcs)
Sanplatec PE Wide Mouth Bottle			02081c	30mL (400pcs)
Sanplatec PE Wide Mouth Bottle			02082c	50mL (400pcs)
Sanplatec PE Wide Mouth Bottle			02083c	100mL (200pcs)
Sanplatec PE Wide Mouth Bottle			02084c	250mL (100pcs)
Sanplatec PE Wide Mouth Bottle			02085c	500mL (100pcs)
Sanplatec PE Wide Mouth Bottle			02086c	1L (50pcs)
Sanplatec PE Wide Mouth Bottle			02087c	2L (30pcs)
Sanplatec PE Wide Mouth Bottle			02088c	3L (30pcs)
Sanplatec PE Wide Mouth Bottle			02089c	5L (15pcs)
Sanplatec PE Wide Mouth Bottle			02090	10L
Sanplatec PE Wide Mouth Bottle			02090c	10L (8pcs)
Sanplatec PE Wide Mouth Bottle (without inner seal)	<ul style="list-style-type: none"> ●Special seal structure ensures ●Cap is molded with transparent ●Smart shape allows easy-grip ●The inner ring of cap fits tight ●Cap is the Polypropylene material ●Main body: PE, Cap: PP ●Without inner seal 	<ul style="list-style-type: none"> ●Each size is for reference only ●The total length includes the 	25010c	100mL (200pcs)
Sanplatec PE Wide Mouth Bottle (without inner seal)			25011c	250mL (100pcs)
Sanplatec PE Wide Mouth Bottle (without inner seal)			25012c	500mL (100pcs)
Sanplatec PE Wide Mouth Bottle (without inner seal)			25013c	1L(50pcs)
Sanplatec PE Wide Mouth Bottle (without inner seal)			25014c	2L(30pcs)
Sanplatec PE Wide Mouth Bottle	<ul style="list-style-type: none"> ●There is dent at the top of b 	<ul style="list-style-type: none"> ●Each size is for reference only ●The height is measured with v 	02091	20L
Sanplatec PE Wide Mouth Bottle (with spigot)	<ul style="list-style-type: none"> ●Spigot thread size: 20A ●With inner seal ●Spigot is made of PP 		02095	5L
Sanplatec PE Wide Mouth Bottle (with spigot)			02096	10L
Sanplatec PE Wide Mouth Bottle (with spigot)			02097	20L
Sanplatec PE Wide Mouth Bottle (with spigot)			02918	5L
Sanplatec PE Wide Mouth Bottle (with spigot)			02919	10L
Square Wash Bottle	<ul style="list-style-type: none"> ●Made in polyethylene ●Efficient square shape to save ●Easy liquid loading with wide ●Tip nozzle is replaceable (so 		02115	250mL
Square Wash Bottle			02116	500mL
Sanplatec Wash Bottle	<ul style="list-style-type: none"> ●Body: LDPE / Cap: HDPE 		02117	250mL (wide Mouth)
Sanplatec Wash Bottle			02118	500mL (Wide Mouth)
Sanplatec Wash Bottle			02119	1L (Wide Mouth)

Book Bottle (Square Bottle) with Spigot	<ul style="list-style-type: none"> ●Book type is efficiently used in a limited space. ●5L,10L: With inner seal (PE) ●Bore diameter: 43mm ϕ ●Spigot thread size (mm): 20A ●Tip/Outer diameter: 17mm ϕ , Inner diameter: 10mm ϕ 		02122	5L PP
Book Bottle (Square Bottle) with Spigot			02123	10L PP
Square Bottle A type	<ul style="list-style-type: none"> ●With inner seal ●Material: Body: PP, Cap: HF 	●Each size is for reference only	02125c	50mL(400pcs)
Square Bottle A type			02126c	100mL (200pcs)
Square Bottle A type			02127c	250mL(100)
Square Bottle A type			02128c	500mL (100pcs)
Square Bottle A type			02129c	1L(50pcs)
Square Bottle B type	<ul style="list-style-type: none"> ●With inner seal ●Body, Inner seal: PE, Cap: P 	●Each size is for reference only	02130c	500mL (100pcs)
Square Bottle B type			02131c	1L(100pcs)
Square Bottle B type			02132c	2L(50pcs)
Square Bottle C type		●Each size is for reference only	02133c	1L(100pcs)
Square Bottle C type		●Gasket is included (not inner)	02134c	2L(50pcs)
Square Bottle (nozzle type)	<ul style="list-style-type: none"> ●Accessory: Plug cap ●Bellows type is expandable and shape memory when bending with being stretched. ●Material: PE 	●Red color has an air ventilation part.	02135	White 2L
Square Bottle (nozzle type)			02136	White 4L
Square Bottle (nozzle type)			02137	Red 2L
Square Bottle (nozzle type)			02138	Red 4L
Sanplatec Rectangular Carboy	●With inner seal		02143	20L
Sanplatec Rectangular Carboy			02145	4L
Sanplatec Rectangular Carboy			02146	5L
Sanplatec Rectangular Carboy			02606	10L Narrow Mouth
Sanplatec Rectangular Carboy			02639	20L Wide Mouth
Sanplatec Rectangular Carboy			02640	20L Narrow Mouth
Sanplatec Rectangular Carboy			26549	30L wide mouth
Sanplatec Sealed Container (round cap)	<ul style="list-style-type: none"> ●Convenient to store and transport ●With inner seal 		02150c	No.1(300pcs)
Sanplatec Sealed Container (round cap)			02151c	No.2(240pcs)
Sanplatec Sealed Container (round cap)			02152c	No.3(120pcs)
Sanplatec Sealed Container (round cap)			02153c	No.4(120pcs)
Sanplatec Sealed Container (round cap)			02154c	No.5(80pcs)
Sanplatec Sealed Container (round cap)			02155c	No.6(80pcs)
Sanplatec Sealed Container (round cap)			02157c	No.8(60pcs)
Sanplatec Sealed Container (round cap)			02158c	No.9(40pcs)
Sanplatec Sealed Container (round cap)	<ul style="list-style-type: none"> ●Convenient to store and transport ●Easy handling for high viscosity ●Lid is flat and stackable. ●Main body/Cap: PP, Inner seal: PE 		02160c	1L(60pcs)
Sanplatec Sealed Container (round cap)			02161c	1.5L(40pcs)
Styrol Round Bottle A type (push lid)		●Be careful with handling because	02186c	5mL(100pcs)
Styrol Round Bottle A type (push lid)			02187c	10mL(100pcs)
Styrol Round Bottle A type (push lid)			02188c	15mL(100pcs)
Styrol Round Bottle A type (push lid)			02189c	25mL(100pcs)

Styrol Round Bottle A type (push lid)			02190c	30mL(100pcs)
Styrol Round Bottle A type (push lid)			02191c	40mL(100pcs)
Styrol Round Bottle A type (push lid)			02192c	50mL(100pcs)
Styrol Round Bottle A type (push lid)			02193c	70mL(50pcs)
Styrol Round Bottle A type (push lid)			02194c	90mL(50pcs)
Styrol Round Bottle A type (push lid)			02195c	120mL(50pcs)
Styrol Round Bottle A type (push lid)			02196c	200mL(50pcs)
Styrol Round Bottle A type (push lid)			02197	300mL
Styrol Round Bottle A type (push lid)			02198	500mL
Styrol Round Bottle C type (screw lid)	●Convenient to store big-size	●Be careful with handling styro	02210c	140mL (50pcs)
Styrol Round Bottle C type (screw lid)			02211c	200mL (50pcs)
Styrol Round Bottle C type (screw lid)			02212	350mL
Styrol Round Bottle C type (screw lid)			02213	600mL
Styrol Square Box	●Applications: Perfect for tra	●Each size is for reference only	02261	No.1(10pcs)
Styrol Square Box			02261c	No.1(50pcs)
Styrol Square Box			02262	No.2(10pcs)
Styrol Square Box			02262c	No.2(50pcs)
Styrol Square Box			02263	No.3(10pcs)
Styrol Square Box			02263c	No.3(50pcs)
Styrol Square Box			02264	No.4(10pcs)
Styrol Square Box			02264c	No.4(50pcs)
Styrol Square Box			02265	No.5(10pcs)
Styrol Square Box			02265c	No.5(50pcs)
Styrol Square Box			02266	No.6(10pcs)
Styrol Square Box			02266c	No.6(55pcs)
Styrol Square Box			02267	No.7(10pcs)
Styrol Square Box			02267c	No.7(35pcs)
Styrol Square Box			02268	No.8(10pcs)
Styrol Square Box			02268c	No.8(28pcs)
Styrol Square Box			02270	No.10(10pcs)
Styrol Square Box			02270c	No.10(50pcs)
Styrol Square Box	●Applications: Perfect for tra	●Each size is for reference only	02271	No.11(10pcs)
Styrol Square Box			02271c	No.11(50pcs)
Styrol Square Box			02272	No.12(10pcs)
Styrol Square Box			02272c	No.12(50pcs)
Styrol Square Box			02273	No.13(10pcs)
Styrol Square Box			02273c	No.13(50pcs)
Styrol Square Box			02274	No.14(10pcs)
Styrol Square Box			02274c	No.14(50pcs)
Styrol Square Box			02275	No.15(10pcs)
Styrol Square Box			02275c	No.15(30pcs)
Styrol Square Box			02276	No.16(10pcs)
Styrol Square Box			02276c	No.16(50pcs)
Styrol Square Box			02278	No.18(10pcs)
Styrol Square Box			02278c	No.18(50pcs)
Styrol Square Box			02279	No.19(10pcs)

Styrol Square Box			02279c	No.19(50pcs)
Styrol Square Box			02280	No.20(10pcs)
Styrol Square Box			02280c	No.20(50pcs)
Styrol Square Box	●Applications: Perfect for tra	●Each size is for reference only	02282	No.22(10pcs)
Styrol Square Box			02282c	No.22(50pcs)
Styrol Square Box			02283	No.23(10pcs)
Styrol Square Box			02283c	No.23(50pcs)
Styrol Square Box			02284	No.24(10pcs)
Styrol Square Box			02284c	No.24(50pcs)
Styrol Square Box			02285	No.25(10pcs)
Styrol Square Box			02285c	No.25(50pcs)
Styrol Square Box			02286	No.26(10pcs)
Styrol Square Box			02286c	No.26(50pcs)
Styrol Square Box			02287	No.27(10pcs)
Styrol Square Box			02287c	No.27(50pcs)
Styrol Square Box			02288	No.28(10pcs)
Styrol Square Box			02288c	No.28(50pcs)
Styrol Square Box			02298	No.21(10pcs)
Styrol Square Box			02298c	No.21(50pcs)
Styrol Square Box			02299	No.17(10pcs)
Styrol Square Box			02299c	No.17(50pcs)
Book Bottle (Square Bottle)	●Book type is efficiently used ●500mL,1L,2L: No inner seal ●5L,10L: With inner seal (PE)	●Caps of 5L and 10L are not as ●5L·10L Cap: Autoclavable	02120	5L PP
Book Bottle (Square Bottle)			02121	10L PP
Book Bottle (square bottle)			02208	1L PP (50pcs)
Book Bottle (square bottle)			02300	500mL PP
Book Bottle (square bottle)			02303	2L PP
Sanplatec RB Bottle	●Hydrocarbons perform excell ●Sanplatec RB bottles are als	●Depending on the liquid used, We recommend to test with a s	02306c	4L Square type (16pcs)
Sanplatec RB Bottle			02307c	2L Square type (40pcs)
Sanplatec RB Bottle			02308c	1L Square type (65pcs)
Sanplatec RB Bottle			02309c	500mL (100pcs)
Sanplatec RB Bottle			02310c	250mL (180pcs)
Sanplatec RB Bottle			02311c	100mL (350pcs)
Square Sample Box (acrylic)	●For storage of important sam ●Square type enables to store		02312	5L
Square Sample Box (acrylic)			02313	10L
Square Sample Box (acrylic)			02314	15L
Square Sample Box (acrylic)			02315	20L
Square Sample Box (acrylic)			02316	30L
Square Sample Box (acrylic)			02317	60L
Styrol Large Square Box	●Suitable for storage of samp ●Packed case dimension: W×	●Styrol is easy to scratch. Plea	02320	O-1
Styrol Large Square Box			02321	O-2
Styrol Large Square Box			02322	O-3
Styrol Large Square Box			02323	O-4
Styrol Large Square Box			02324	O-5
Styrol Large Square Box			02325	O-6
Styrol Large Square Box			02326	O-7

Styrol Large Square Box		02327	O-8
Round Sample Box (acrylic)	<ul style="list-style-type: none"> ● Fully transparent and convenient ● Easy view from any angle. 	02340	5L
Round Sample Box (acrylic)		02341	10L
Round Sample Box (acrylic)		02342	15L
Round Sample Box (acrylic)		02343	20L
Round Sample Box (acrylic)		02344	30L
Round Sample Box (acrylic)		02345	60L
TPX Dish	<ul style="list-style-type: none"> ● Petri dishes molded with TPX ● Autoclavable :121°C ● 2 size type: shallow, deep 	02360c	Shallow type (10pcs)
TPX Dish		02361c	Deep type (10pcs)
Sanplatec Square Tray	<ul style="list-style-type: none"> ● The bottom is flat. 	02370	S-1
Sanplatec Square Tray		02371	S-2
Sanplatec Square Tray		02372	S-3
Large Size Tray	<ul style="list-style-type: none"> ● The bottom surface is irregular 	02373	Type A
Large Size Tray		02374	Type B
Large Size Tray		02375	Type C
Large Size Tray		02376	Type D
Large Size Tray		02377	Type E
Large Size Tray		02378	Type F
Large Size Tray		02379	Type G
Bottle Tray	<ul style="list-style-type: none"> ● Tray: polyethylene /375 × 275 ● Holder: vinyl chloride ● For 6pcs of 500mL bottle / 12 	02387	2 type
Bottle Tray		02387c	2 type(5pcs)
Bottle Tray	<ul style="list-style-type: none"> ● Tray: polyethylene/375 × 275 ● Holder: vinyl chloride ● For 24pcs of 25mL bottle / 12 	02388	3 type
Bottle Tray		02388c	3 type(5pcs)
Square Tray	<ul style="list-style-type: none"> ● Outer dimension (mm): 390 × 390 ● Material/ vinyl chloride 	02389	1 type
Square Tray		02389c	1 type(5pcs)
Square Tray	<ul style="list-style-type: none"> ● Outer dimension(mm):440 × 440 ● Material/ vinyl chloride 	02390	1 type
Square Tray		02390c	1type (5pcs)
Tabletop Bottle Holder	<ul style="list-style-type: none"> ● Hole diameter: 85mm ● Spec: 6pcs of 500mL bottle ● Outer dimension of main body: 180 × 180 × 100 	02391	TH-6
Tabletop Bottle Holder	<ul style="list-style-type: none"> ● Cylinder I.D. 82mm × 80H (All type) 	02393	2 type
Tabletop Bottle Holder		02394	3 type
Tabletop Bottle Holder		02395	4 type
Tabletop Bottle Holder		02396	6 type
Bottle Protector (Safety Netron)	<ul style="list-style-type: none"> ● Perfect for protection of reaction vessels ● Stretch material can fit any size ● Excellent cold resistance, no deformation ● Easy cutting with scissors and utility knife ● Scale (m): 5 ● Material: Ethylene-vinyl acetate 	02400	N-1 (5M)
Bottle Protector (Safety Netron)		02401	N-2 (5M)
Bottle Protector (Safety Netron)		02402	N-3 (5M)
Bottle Protector (Safety Netron)		02403	N-4 (5M)
Bottle Protector (Safety Netron)		02404	N-5 (5M)
Bottle Protector (Safety Netron)		02405	N-6 (5M)
Styrol Round Box	<ul style="list-style-type: none"> ● There is a transparent cover ● The lid contains small holes ● Stackable with filling water ● Material: Aquatic tanks PS/PE 	02441	R1 type φ 240MM

Styrol Round Box			02442	R2 type φ 300MM
Styrol Round Box			02443	R1 type Lid
Styrol Round Box			02444	R2 type Lid
Chemical Storage Shelf	<ul style="list-style-type: none"> ●Main body: 10mm low-foam ●In both the upper and lower ●The upper row has a key and ●The bottom row has drawers ●A key is also attached to the 	<ul style="list-style-type: none"> ●There is no resistance to corr 	02454	P600 type (Upper part)
Chemical Storage Shelf			02455	P600 type (Lower part)
Chemical Storage Shelf			02456	P600 type set
Chemical Storage Shelf	<ul style="list-style-type: none"> ●Main body: 10mm low-foam ●In both the upper and lower ●The upper row has a key and ●The bottom row has drawers ●The door at the bottom row 	<ul style="list-style-type: none"> ●There is no resistance to corr ●Withstand weight of shelves a 	02457	P900 type (Upper part)
Chemical Storage Shelf			02458	P900 type (Lower part)
Chemical Storage Shelf			02459	P900 type set
Sanplatec PE Narrow Mouth Bottle (with handle)	<ul style="list-style-type: none"> ●With inner seal ●Handle material:PP 		02505	5L
Sanplatec PE Narrow Mouth Bottle (with handle)			02506	10L
Sanplatec PE Wide Mouth Bottle (with handle)	<ul style="list-style-type: none"> ●With inner seal ●Handle material:PP 		02507	5L
Sanplatec PE Wide Mouth Bottle (with handle)			02508	10L
Sanplatec PE Wide Mouth Bottle (with handle)			02924	5L
Sanplatec PE Wide Mouth Bottle (with handle)			02925	10L
Chemical Storage Cabinet	<ul style="list-style-type: none"> ●Convenient for the usage on ●Easy view through light-resis ●The bottom part is bat shape ●Lockable for safety ●Possible to store max 15pcs 		02535	450 Atype (Upper part)
Chemical Storage Cabinet	<ul style="list-style-type: none"> ●Convenient for usage on a lab table. ●The front door can be also used as a worktable. ●The bottom part is bat shape. Safe when a small amount of chemical leakage occurs. ●Lockable for safety ●Possible to store max 15pcs of 500mL bottle 		02536	450 Btype (Lower part)
PTFE Tube Rack	<ul style="list-style-type: none"> ●Assemble type. Easy for clea ●Excellent heat and cold resis ●Can be used in a freezer or ●Autoclavable 		02560	TR-1
PTFE Tube Rack			02561	TR-3
PTFE Tube Rack			02562	TR-2
PTFE Tube Rack			02563	TR-4
Mouse/Rat Restrainer (capsule type)	<ul style="list-style-type: none"> ●This is made of transparent ●Partition plate allows the sp 	<ul style="list-style-type: none"> ●When cleaning, wipe off with a ●Cleaning and disinfection with It may cause material deteriora 	02564	RC-1 type
Mouse/Rat Restrainer (capsule type)			02565	RC-2 type
Mouse/Rat Restrainer (capsule type)			02566	RC-3 type
Animal Collar (mouse type)	<ul style="list-style-type: none"> ●Easy close with velcro strap ●Protection of the application ●Perfect fit for animal neck a 		02580	(10pcs)
Sanplatec Rectangular Carboy with Spigot	<ul style="list-style-type: none"> ●PP spigot: 25A ●With inner seal 		02641	20L Narrow Mouth
Sanplatec Rectangular Carboy with Spigot			02642	20L Wide Mouth
PE Wide Mouth Amber Bottle (with inner seal)	<ul style="list-style-type: none"> ●Excellent light-shielding perf ●Ideal for storing ultraviolet re ●5L, 10L has a hole through v 	<ul style="list-style-type: none"> ●Each size is for reference only ●Total length includes a cap. ●With inner seal (white) 	02910c	100mL (200pcs)
PE Wide Mouth Amber Bottle (with inner seal)			02911c	250mL (100pcs)
PE Wide Mouth Amber Bottle (with inner seal)			02912c	500mL (100pcs)

PE Wide Mouth Amber Bottle (with inner seal)		02913c	1L (50pcs)
PE Wide Mouth Amber Bottle (with inner seal)		02914c	2L (30pcs)
PE Wide Mouth Amber Bottle (with inner seal)		02915c	3L (30pcs)
PE Wide Mouth Amber Bottle (with inner seal)		02916c	5L (15pcs)
PE Wide Mouth Amber Bottle (with inner seal)		02917c	10L (8pcs)
PE Precipitation Bottle	<ul style="list-style-type: none"> ●With inner seal ●Position of stopcock is 1/3 ●Stopcock: PP 	02995	5L Wide Mouth
PE Precipitation Bottle		02996	10L Wide Mouth
PE Precipitation Bottle		02997	20L Wide Mouth
PP Spigot	<ul style="list-style-type: none"> ●Dedicated spigot for Sanplat ●Rubber seal is viton rubber. ●Please use a hose with inner 	03040	20A
PP Spigot		03041	25A
PE Spigot	<ul style="list-style-type: none"> ●Dedicated spigot for Sanplat ●Rubber seal is viton rubber. 	03043	20A
PE Spigot		03044	25A
Labmate Arm Cover	<ul style="list-style-type: none"> ●Whole length : 370mm ●U-1 type includes 6 cases of 	03050	U-1(300pcs)
Labmate Arm Cover		03051	U-2(50pcs)
Dispensing Bottle (spray type)	●Bottle material:PET	03489	S-50
Dispensing Bottle (spray type)		03490	S-150
Dispensing Bottle (spray type)		03491	S-300
Dispensing Bottle (spray type)		03492	S-500
Dispensing Bottle (pump type)	●1 pump: approx. 1mL	03493	P-150
Dispensing Bottle (pump type)		03494	P-300
Dispensing Bottle (pump type)		03495	P-500
PP Tube Connector (taper type)		03530	Itype Large (10pcs)
PP Tube Connector (taper type)		03531	Itype Middle (10pcs)
PP Tube Connector (taper type)		03532	Itype Small (10pcs)
PP Tube Connector (taper type)		03533	Different diameter type Large (10pcs)
PP Tube Connector (taper type)		03534	Different diameter type Small (10pcs)
PP Tube Connector (taper type)		03535	Y type Large (10pcs)
PP Tube Connector (taper type)		03536	Y type Middle (10pcs)
PP Tube Connector (taper type)		03537	Y type Small (10pcs)
PP Tube Connector (taper type)		03538	Ttype Large (10pcs)
PP Tube Connector (taper type)		03539	Ttype Middle (10pcs)
PP Tube Connector (taper type)		03540	Ttype Small (10pcs)
TPX® Tube Connector (taper type)	<ul style="list-style-type: none"> ●Great chemical and heat resistance, and durability. ●Heat resistance 170°C. Autoclave sterilization is possible. ●Transparent enables easy viewing of fluid. 	03541	Itype Large (10pcs)
TPX® Tube Connector (taper type)		03542	Itype Middle (10pcs)
TPX® Tube Connector (taper type)		03543	Itype Small (10pcs)
TPX® Tube Connector (taper type)	<ul style="list-style-type: none"> ●Great chemical and heat resistance, and durability. ●Heat resistance 170°C. Autoclave sterilization is possible. ●Transparent enables easy viewing of fluid. 	03544	Different diameter Large (10pcs)
TPX® Tube Connector (taper type)		03545	Different diameter Small (10pcs)

TPX® Tube Connector (taper type)	<ul style="list-style-type: none"> ● Great chemical and heat res ● Heat resistance 170°C. Aut ● Transparent enables easy vi 		03546	Y type Large (10pcs)
TPX® Tube Connector (taper type)			03547	Y type Middle (10pcs)
TPX® Tube Connector (taper type)			03548	Y type Small (10pcs)
TPX® Tube Connector (taper type)	<ul style="list-style-type: none"> ● Great chemical and heat res ● Heat resistance 170°C. Aut ● Transparent enables easy vi 		03549	Ttype Large (10pcs)
TPX® Tube Connector (taper type)			03550	Ttype Middle (10pcs)
TPX® Tube Connector (taper type)			03551	Ttype Small (10pcs)
Vacuum Filter	<ul style="list-style-type: none"> ● Example of application: radio ● Outer dimension (mm):235 x ● Filter size: 25mm φ ● Please use silicon plug No. 5 ● The number of processing se ● Well bore diameter: 18mm φ ● Hose nipple outer diameter: ● O-ring: silicon rubber ● VF-10N enables to collect fi ● Easy and simple operation. F 	● Filters are glass fiber filter pa	04000	VF-10N
Vacuum filter			04034	VF-10
Sample Cup	<ul style="list-style-type: none"> ● Select according to material 		04036	1A (2000pcs)
Sample Cup			04037	2A (2000pcs)
Sample Cup			04038	3A (1000pcs)
Sample Cup			04040	5A (1000pcs)
Sample Cup			04041	5A-1 (1000pcs)
Sample Cup			04042	6A (1000pcs)
Sample Cup			04043	7A (1000pcs)
Sample Cup			04044	8A (1000pcs)
Sample Cup			04045	9A (1000pcs)
Sample Cup			04048	12A (1000pcs)
Sample Cup			04049	13A (2000pcs)
Test Tube	<ul style="list-style-type: none"> ● Sterlized item: EOG sterilizat 		04121	9-D (1000pcs)
Test Tube			04122	10-D (1000pcs)
Test Tube			04123	11-D (1000pcs)
Test Tube			04125	13-D (1000pcs)
Test Tube			04126	14-D (1000pcs)
Test Tube			04127	15-D (1000pcs)
Test Tube			04128	16-D (1000pcs)
Test Tube			04129	17-D (1000pcs)
Test Tube			04130	9-D (1000pcs)
Test Tube			04131	10-D (1000pcs)
Test Tube			04132	11-D (1000pcs)
Test Tube			04134	13-D (1000pcs)
Test Tube			04135	14-D (1000pcs)
Test Tube			04136	15-D (1000pcs)
Test Tube			04137	16-D (1000pcs)
Test Tube			04138	17-D (1000pcs)
PTFE Connector (T type)			04170	PT-1
PTFE Connector (T type)			04171	PT-2
PTFE Connector (T type)			04172	PT-3
PTFE Connector (T type)			04173	PT-4

PTFE Connector (L type)			04174	PE-1 (1/4)
PTFE Connector (L type)			04175	PE-2 (3/8)
PTFE Connector (L type)			04176	PE-3 (1/2)
PTFE Connector (L type)			04177	PE-4
PTFE Connector (nipple type)			04182	N-1
PTFE Connector (nipple type)			04183	N-2
PTFE Connector (nipple type)			04184	N-3
PTFE Connector (nipple type)			04185	N-4
PTFE Connector (nipple type)			04186	N-5
PTFE Connector (nipple type)			04187	N-6
PTFE Connector (nipple type)			04188	N-7
PTFE Connector (nipple type)			04189	N-8
PTFE Connector (nipple type)			04190	N-9
PTFE Connector (tube/nipple type)			04191	HM-1
PTFE Connector (tube/nipple type)			04192	HM-2
PTFE Connector (tube/nipple type)			04193	HM-3
PTFE Connector (tube/nipple type)			04194	HM-4
PTFE Connector (tube/nipple type)			04195	HM-5
PTFE Connector (tube/nipple type)			04196	HF-1
PTFE Connector (tube/nipple type)			04197	HF-2
PTFE Connector (tube/nipple type)			04198	HF-3
PTFE Connector (tube/nipple type)			04199	HF-4
PTFE Connector (tube/nipple type)			04200	HF-5
PTFE Connector (coupling type)			04201	CL-1
PTFE Connector (coupling type)			04202	CL-2
PTFE Connector (coupling type)			04204	MF-1
Sanplatec PTFE Stopper			04209	SP-1
Sanplatec PTFE Stopper			04210	SP-2
Sanplatec PTFE Stopper			04211	SP-3
Sanplatec PTFE Stopper			04212	SP-4
Sanplatec PTFE Stopper			04213	SP-5
Sanplatec PTFE Stopper			04214	SP-6
PTFE Mortar (with pestle)	<ul style="list-style-type: none"> ● With PTFE pestle ● Thickness 5mm (side) 		04234	DD-1(PTFE)
PTFE Mortar (with pestle)			04235	DD-2(PTFE)
PTFE Mortar (with pestle)			04236	DD-3(PTFE)
PTFE Mortar (with pestle)			04237	DD-4(PTFE)
PTFE Mortar (with pestle)			04238	DD-5(PTFE)
Protector with Gloves	<ul style="list-style-type: none"> ● Made of acrylic, Excellent transparent. ● Convenient for protection when testing ● Acrylic board: thickness 5mm 	<ul style="list-style-type: none"> ● Make sure to fix on your worktable when using. ● Not for RI protection. 	04307	GP-1
Protector with Gloves			04308	GP-2
Protector with Gloves			04309	GP-3
PFA Square Container (with lid)	<ul style="list-style-type: none"> ● This is manufactured by inte ● Use as a bat for reagent exp ● Adaptable slide size: 25x75m 	<ul style="list-style-type: none"> ● Container depth is shallow. Pl 	04420	38X82X19

ETFE Narrow Mouth Bottle	<ul style="list-style-type: none"> ●Cap: Polyethylene ●Thickness: 2mm (for reference) ●ETFE has the similar characteristic as PTFE. There is some resistance for acid, acids, alkalis, and organic solvents. ●Heat resistant temperature approx. 150°C ●Excellent mechanical strength, electrical insulation, and radiation resistance 		04922	5L
ETFE Narrow Mouth Bottle			04923	10L
ETFE Narrow Mouth Bottle			04924	20L
ETFE Narrow Mouth Bottle			04925	30L
ETFE Wide Mouth Bottle	<ul style="list-style-type: none"> ●Cap: Polyethylene ●Thickness(mm):2 (reference) ●ETFE has the similar characteristic as PTFE. ●Heat resistant temperature approx. 150°C ●Excellent mechanical strength, electrical insulation, and radiation resistance 		04926	5L
ETFE Wide Mouth Bottle			04927	10L
ETFE Wide Mouth Bottle			04928	20L
ETFE Wide Mouth Bottle			04929	30L
PTFE Test Tube (round bottom)			04934	8mL
PTFE Test Tube (round bottom)			04935	17mL
PTFE Test Tube (round bottom)			04936	8mL (with cap)
PTFE Test Tube (round bottom)			04937	17mL (with cap)
PTFE Test Tube (conical type)	<ul style="list-style-type: none"> ●Test tube made of fluorinated resin ●Convenient to concentrate sample at the bottom. ●Two types are available: with cap and without cap ● Excellent chemical resistance ●Not transparent. No clear viewing of contents. 		04938	13mL (with cap)
PTFE Test Tube (conical type)			04939	15mL
PTFE Test Tube (conical type)			04940	13mL (with cap)
PTFE Test Tube (conical type)			04941	15mL (with cap)
PTFE Tweezers	<ul style="list-style-type: none"> ●A guide rod is attached so that the sample can be easily taken out. 		04957	100L (with guide)
PTFE Tweezers			04958	150L (with guide pin)
PTFE Tweezers			04959	200L (with guide pin)
PTFE Tweezers	<ul style="list-style-type: none"> ●Tip shape: flat 		04960	100L (with guide pin)
PTFE Tweezers	<ul style="list-style-type: none"> ●A guide rod is attached so that the sample can be easily taken out. 		04961	150L (with guide pin)
PTFE Tweezers			04962	200L (with guide pin)
PTFE Stopper	<ul style="list-style-type: none"> ●Great chemical resistance ●Excellent sealing performance 		04967	14/23
PTFE Stopper			04968	19/26
PTFE Stopper			04969	24/29
PTFE Stopper			04970	29/32
PTFE Stopper			04971	34/35
PFA Erlenmeyer Flask	<ul style="list-style-type: none"> ●High temperature resistance ●Flat bottom supports stability 		05027	PFA 100mL
PFA Erlenmeyer Flask			05028	PFA 200mL
PFA Dropping Pipette			05034	2.5mL
PFA Dropping Pipette			05035	5mL

Gas Sampling Bag (Tedlar® Bag)	<ul style="list-style-type: none"> ● Suction port is outer diameter ● The Tedlar® bag uses Dupont ● Tedlar® is a registered brand ● Excellent gas barriers can effectively ● It does not contain any plasticizers ● It is used for air pollution control ● Also used for sampling for air 	05311	1L
Gas Sampling Bag (Tedlar® Bag)		05312	2L
Gas Sampling Bag (Tedlar® Bag)		05313	3L
Gas Sampling Bag (Tedlar® Bag)		05314	5L
Gas Sampling Bag (Tedlar® Bag)		05315	10L
Gas Sampling Bag (Tedlar® Bag)		05316	15L
Gas Sampling Bag (Tedlar® Bag)		05317	20L
Gas Sampling Bag (Tedlar® Bag)		05318	30L
Gas Sampling Bag (Tedlar® Bag)		05319	50L
Gas Sampling Bag (Tedlar® Bag)		05320	100L
Gas Sampling Bag (Tedlar® Bag)		05321	Stopper TS-10 (10pcs)
PTFE Tube	<ul style="list-style-type: none"> ● Excellent non-adhesiveness ● Suitable for electrical insulation ● Excellent weather tolerance 	05450	21 (10M)
PTFE Tube		05451	3 × 2 (10M)
PTFE Tube		05452	4 × 3 (10M)
PTFE Tube		05453	5 × 4 (10M)
PTFE Tube		05454	6 × 5 (10M)
PTFE Tube		05455	7 × 6 (10M)
PTFE Tube		05456	8 × 7 (10M)
PTFE Tube		05457	9 × 8 (10M)
PTFE Tube		05458	10 × 9 (10M)
PTFE Tube		05459	11 × 10 (10M)
PTFE Tube		05460	12 × 11 (10M)
PTFE Tube		05461	13 × 12 (10M)
PTFE Tube	<ul style="list-style-type: none"> ● Excellent non-adhesiveness ● Suitable for electrical insulation ● Excellent weather tolerance 	05462	3 × 1 (10M)
PTFE Tube		05463	4 × 2 (10M)
PTFE Tube		05464	5 × 3 (10M)
PTFE Tube		05465	6 × 4 (10M)
PTFE Tube		05466	7 × 5 (10M)
PTFE Tube		05467	8 × 6 (10M)
PTFE Tube		05468	9 × 7 (10M)
PTFE Tube		05469	10 × 8 (10M)
PTFE Tube		05470	11 × 9 (10M)
PTFE Tube		05471	12 × 10 (10M)
PTFE Tube		05472	13 × 11 (10M)
PTFE Tube		05473	14 × 12 (10M)
PTFE Tube		05474	15 × 13 (10M)
PTFE Tube		05475	16 × 14 (10M)
PTFE Tube		05478	18 × 16 (10M)
PTFE Tube		05479	21 × 19 (10M)
PFA Tube	<ul style="list-style-type: none"> ● Excellent non-adhesiveness ● Suitable for electrical insulation ● Excellent weather tolerance 	05657	3 × 2
PFA Tube		05658	4 × 2
PFA Tube		05658	4 × 2

PFA Tube			05659	4 × 3
PFA Tube			05660	6 × 4
PFA Tube			05661	6 × 5
PFA Tube			05662	8 × 6
PFA Tube			05663	8 × 7
PFA Tube			05664	10 × 8
PFA Tube			05666	12 × 10
PFA Tube			05668	14 × 12
PFA Tube			05669	15 × 13
PFA Tube			05670	16 × 13
PFA Tube			05671	16 × 14
PFA Tube			05672	18 × 15
PFA Tube			05673	18 × 16
PFA Tube			05674	19 × 16
PFA Tube			05675	22 × 19
PFA Tube			05677	25 × 22
PFA Tube	<ul style="list-style-type: none"> ●Excellent non-adhesiveness ●Suitable for electrical insulation ●Excellent weather tolerance 		05679	3 × 2 × 0.5
PFA Tube			05680	4 × 2 × 1.0
PFA Tube			05681	4 × 3 × 0.5
PFA Tube			05682	6 × 4 × 1.0
PFA Tube			05683	6 × 5 × 0.5
PFA Tube			05684	8 × 6 × 1.0
PFA Tube			05685	8 × 7 × 0.5
PFA Tube			05686	10 × 8 × 1.0
PFA Tube			05688	12 × 10 × 1.0
PFA Tube			05690	14 × 12 × 1.0
PFA Tube			18470	15 × 13 × 1.0
PFA Tube			18471	16 × 13 × 1.5
PFA Tube			18472	16 × 14 × 1.0
PFA Tube			18473	18 × 15 × 1.5
PFA Tube			18474	18 × 16 × 1.0
PFA Tube			18475	19 × 16 × 1.5
PFA Tube			18476	22 × 19 × 1.5
PFA Tube			18477	25 × 22 × 1.5
Sanpail	<ul style="list-style-type: none"> ●Plastic pail can ●Convenient for temporary storage ●Take a lid off easily using a screwdriver 	<ul style="list-style-type: none"> ●Body and lid are separately stored 	05819	NO.6 Main Body
Sanpail			05820	NO.6 Lid
Sanpail			05821	NO.12 Main Body
Sanpail			05822	NO.12 Lid
Sanpail			05823	NO.20 Main Body
Sanpail			05824	NO.20 Lid
Glove Box (gas purge type)	<ul style="list-style-type: none"> ●All acrylic material enables to perform gas analysis ●Body and lid are detachable ●Pass box can perform gas exchange ●Accessory: vacuum gauge for pressure measurement ●Pressure resistant: Main body 	<ul style="list-style-type: none"> ●Outer dimension of main body 	06103	DOM-1
Glove Box (gas purge type)			06104	DOM-2

Sink (PVC two tanks)	<ul style="list-style-type: none"> ●Excellent acid and alkali tolerance ●The structure of inner bottom allows less fluid remaining in sink. ●Drain outlet: 25A hose nipple ●Comes with a drain hose 3m 		06501	N900
Sink (PVC Two Tanks)			06502	N1200
Cylinder (small scale)	<ul style="list-style-type: none"> ●TPX® is a registered™ of Mits 		06503	20mL(PP)
Cylinder (small scale)			06504	100mL(PP)
Cylinder (small scale)			06505	20mL(TPX)
Cylinder (small scale)			06506	100mL(TPX)
Sink (PVC one tank)	<ul style="list-style-type: none"> ●Excellent acid and alkali tolerance ●The structure of inner bottom allows less fluid remaining in sink. ●Sink depth: 150mmH ●Drain outlet: 25A hose nipple ●Comes with a drain hose 3m 		06507	600 type
Sink (PVC one tank)			06508	900 type
Sink (PVC one tank)			06509	1200 type
Graduated Beaker with Handle	<ul style="list-style-type: none"> ●Great transparent ●Heat resistance temperature ●Agitation with magnetic stirrer ●If chemical resistance and he 		06510	200mL
Graduated Beaker with Handle			06511	300mL
Graduated Beaker with Handle			06512	500mL
Graduated Beaker with Handle			06513	1L
Graduated Beaker with Handle			06514	2L
Sanplatec Wash Bottle	<ul style="list-style-type: none"> ●Nozzle and cap are integrate 		06515	125mL
Sanplatec Wash Bottle			06516	250mL
Sanplatec Wash Bottle			06517	500mL
Sanplatec Wash Bottle			06518	1L
Sanplatec PE Wide Mouth Bottle (with spigot)	<ul style="list-style-type: none"> ●Spigot thread size:20A ●With inner seal ●Handle and spigot : PP 		02926	5L
Sanplatec PE Wide Mouth Bottle (with spigot)			02927	10L
Sanplatec PE Wide Mouth Bottle (with spigot)			06521	5L
Sanplatec PE Wide Mouth Bottle (with spigot)			06522	10L
Acrylic Vacuum Glove Box Clear View	<ul style="list-style-type: none"> ●Easy handling to move items ●Use the right pass box to take ●Accessory Valve: body 2pcs, pass box 2 Vacuum gauge: body+pass be Glove+Glove band: each 1 pair Glove port door: 1pair (when ●Valve nozzle diameter: φ9mm ●Internal dimension of pass box 	<ul style="list-style-type: none"> ●Do not operate with sudden c ●Vacuum gauge shows the pres ●Do not operate when the insid ●The outer dimension of main b 	06590	VG-800 type
Glove Box (Gas Purge type)	<ul style="list-style-type: none"> ●Excellent clear view ●Pass box corresponds to the ●Accessory: glove 1pair, need ●Open size of door attached i 	<ul style="list-style-type: none"> ●Main body part cannot be dec ●The outer dimension of main b 	06591	DX-II type
Draft Tabletop Type	<ul style="list-style-type: none"> ●Made of vinyl chloride with e ●Just place on lab table or we ●Duct :100mm φ x3m, Water s < Spec of fan > ●100V, 50/60(Hz) ●Input(W): 16/15 ●Rotation frequency(R.P.M): 2 ●Max air volume(m³/min): 2.7 ●Max static pressure(mm/H₂O) ●Noise(db(A)): 42/46 	<ul style="list-style-type: none"> ●Do not use flammable substar ●Fans and fittings are not chim 	06592	DF-1
Acrylic Vacuum Glove Box Clear View	<ul style="list-style-type: none"> ●Spec: main body: acrylic(tra ●Pass box: acrylic(transparent ●Vacuum valve: R1/4 nozzle: 	<ul style="list-style-type: none"> ●Do not operate when the insid ●The outer dimension of main b 	06593	M-600H type

Sanplatec Pre-Cleaned Wide Mouth PE Bottle	●We clean Sanplatec PE Wide Mouth Bottle with pure water in clean rooms (cleanliness: class 100)and double packed in class 1000 environment.		06676	250mL(5pcs)
Sanplatec Pre-Cleaned Wide Mouth PE Bottle			06677	100mL(5pcs)
Sanplatec Pre-Cleaned Wide Mouth PE Bottle			06678	500mL(2pcs)
Sanplatec Pre-Cleaned Wide Mouth PE Bottle			06679	1L(1pc)
PTFE Spatula (trapezoid type)	●Excellent heat and chemical		06723	50 × 40mm
PTFE Spatula (trapezoid type)			06724	70 × 50mm
PTFE Spatula (trapezoid type)			06725	90 × 60mm
PTFE Spatula (trapezoid type)			06726	120 × 85mm
PTFE Spatula (semicircle type)	●Excellent heat and chemical		06727	75 × 45mm
PTFE Spatula (semicircle type)			06728	100 × 55mm
PTFE Spatula (semicircle type)			06729	125 × 75mm
PTFE Spatula (semicircle type)			06730	150 × 85mm
PFA Graduated Cylinder	●Excellent heat and chemical resistance ●Cylinder part: PFA ●Base part: PTFE		06766	10mL
PFA Graduated Cylinder			06767	30mL
PFA Graduated Cylinder			06768	50mL
PFA Graduated Cylinder			06769	100mL
PFA Graduated Cylinder			06770	200mL
PFA Graduated Cylinder			06771	300mL
PFA Graduated Cylinder			06772	500mL
PFA Graduated Cylinder			06773	1L
PFA Graduated Cylinder			06774	2L
Glove Box (Gas Purge type)	●G and GS type are available. ●Bottom part of GS type is made of stainless steel. ●Pass box perform to decompress. ●The inner dimension of pass box: ●Stopcock for gas purge mouth. ●Valve on pass box: 2pcs(no. 1000). ●Vacuum gauge on pass box: 1pc. ●Glove and Glove band: 1 set. ●Material of bottom board: acrylic.	●Main body is not allowed to decompress. ●The outer dimension of main body is 250mm × 250mm × 150mm.	06903	G type
Acrylic Vacuum Glove Box Clear View	●VSC-1000 type allows you to see the inside of the box. ●Front opening size: 250mm × 250mm. ●Easy gas purge operation via pass box. ●The inner diameter of pass box: 100mm. ●Valve on main body and pass box: 2pcs. ●Glove and Glove band: 1 set. ●Glove port door: 1 pair (wheel type). <Unit (control box)> Spec ●Vacuum gauge: main body, 1pc ●Pressure gauge: 1pc (0~1MPa) ●Regulator: 1pc ●Flow meter: main body 1pc(200L/min) ●Weight: approx. 6kg	●Vacuum pump is separately set. ●Do not operate when the inside of the box is decompressed. ●The outer dimension of main body is 250mm × 250mm × 150mm.	06904	VSC-1000 type (with Unit)
PFA Dish			08136	50
PFA Dish			08137	62
PFA Dish			08138	76
PFA Dish			08139	100
PFA Funnel	●Outer diameter (mm): 91 ●Total length (mm): 114 ●Stem: 10φ × 15φ × 60L(mm) ●Thickness: 2.5mm		08424	90mm
Sanplatec PFA Test Tube	●Excellent chemical resistance ●Excellent heat and cold resistance ●Non-adhesive. Worry-free with sticky substances. ●Translucent enables easy view		03498	10mL
Sanplatec PFA Test Tube			03499	15mL
Sanplatec PFA Test Tube			08553	5mL

Sanplatec PFA Test Tube			08556	20mL
Sanplatec PFA Test Tube			08557	25mL
Sanplatec PFA Test Tube			08558	50mL
Sanplatec PFA Test Tube			08559	60mL
Sanplatec PFA Graduated Beaker with Handle	<ul style="list-style-type: none"> ●Integrated molding of main body ●Useful in oilbath ●The handle fits comfortably ●The inside bottom part is flat ●The protruding scales assist ●All capacities are marked with ●Excellent weather resistance ●Continuous usage OK 		11110	100mL
Sanplatec PFA Graduated Beaker with Handle			11111	200mL
Sanplatec PFA Graduated Beaker with Handle			11112	500mL
Sanplatec PFA Graduated Beaker with Handle			11113	1L
Sanplatec PFA Graduated Beaker with Handle			11114	2L
Sanplatec PFA Graduated Beaker	<ul style="list-style-type: none"> ●The inside bottom part is flat ●The protruding scales assist ●All capacities are marked with ●PFA material has non-adhesive ●Excellent weather resistance ●Continuous usage OK 		11115	30mL
Sanplatec PFA Graduated Beaker			11116	50mL
Sanplatec PFA Graduated Beaker			11117	100mL
Sanplatec PFA Graduated Beaker			11118	200mL
Sanplatec PFA Graduated Beaker			11119	300mL
Sanplatec PFA Graduated Beaker			11120	500mL
Sanplatec PFA Graduated Beaker			11121	1L
Mouse Restrainer with defecation	<ul style="list-style-type: none"> ●Easy operation with a base t ●It can hold tight mouse's head ●A stopper is attached not to ●Mouse holder size (mm): $\phi 3$ 		13302	Mouse KA
Rat Restrainer with defecation	<ul style="list-style-type: none"> ●Easy operation with a base t ●It can hold tight mouse's head ●A stopper is attached not to ●Mouse holder size (mm): $\phi 7$ 		13303	Rat KK
Polystyrene Dish (electron beam sterilized)	<ul style="list-style-type: none"> ●Size(mm): $\phi 90 \times 15$ ●Electron beam sterilized 		13395	$\phi 90 \times 15$ (10pcs/bag)
Polystyrene Dish (square type)	<ul style="list-style-type: none"> ●Electron beam sterilized 		13868	Square type (240pcs)
PFA Funnel	<ul style="list-style-type: none"> ●Angle($^{\circ}$): 60 ●Stem (outer diameter \times inner) 		15236	5143
PFA Funnel			15237	5146
Sanplatec PFA Spoon			15239	174 \times 22 \times 10H mm
PFA Test Tube (Screw Cap type)	<ul style="list-style-type: none"> ●Excellent chemical resistance ●Heat and cold resistance. ●Translucent color allows you ●Very low eluted metal ion. It ●Non-sticky surface. Suitable 		15398	8215R
PFA Test Tube (Screw Cap type)			15399	8215C
PFA Test Tube (Screw Cap type)			15400	8226R
PFA Test Tube (Screw Cap type)			15401	8232R
PFA Test Tube (Screw Cap type)			15402	8250R
PFA Test Tube (Screw Cap type)			15403	8250C
PFA Test Tube (Screw Cap type)			15404	8260
Sanplatec PFA Square Container	<ul style="list-style-type: none"> ●Large size of PFA square container 	<ul style="list-style-type: none"> ●The lid is separately sold. 	15415	1.6 Main Body
Sanplatec PFA Square Container			15416	4.3 Main Body
Sanplatec PFA Square Container			15417	Lid only

Sanplatec PFA Tray	<ul style="list-style-type: none"> ● Great chemical resistance. Not suitable to soak semiconductor ● Suitable to soak semiconductor ● Available in a wide range of sizes ● Great non-adhesiveness 	● There is unevenness at the inner surface	15418	Small type S
Sanplatec PFA Tray			15419	Small type M
Sanplatec PFA Tray			15420	Small type L
Sanplatec PFA Tray			15421	C-1
Sanplatec PFA Tray			15422	C-2
Sanplatec PFA Tray			15423	C-3
Sanplatec PFA Tray			15424	C-4
Sanplatec PFA Tray			15425	C-5
Sanplatec PFA Tray			15426	BS
Sanplatec PFA Tray			15427	BM
Sanplatec PFA Tray			15428	BL
Sanplatec PFA Tank	<ul style="list-style-type: none"> ● This tank is large and has deep tank ● This can be also used as a container ● You can also use a basket set 		15431	Set (Tank & Basket)
Sanplatec PFA Tank			15432	Tank only
Sanplatec PFA Tank			15433	Basket only
PFA Graduated Pipette	● All part is made of PFA.		16156	2mL
PFA Graduated Pipette			16157	5mL
PFA Graduated Pipette			16158	10mL
PFA Round Bottom Flask	<ul style="list-style-type: none"> ● Central taper size: 29/42 ● Side taper size: 15/35 	● This product can be used in a vacuum chamber. When used in decompression suit	16159	500mL
PFA Round Bottom Flask			16160	500mL
PFA Round Bottom Flask			16161	500mL
PTFE Pipe Bushing			16241	PB-1-2
PTFE Pipe Bushing			16242	PB-1-3
PTFE Pipe Bushing			16243	PB-2-3
PTFE Pipe Bushing			16244	PB-2-4
PTFE Pipe Bushing			16245	PB-3-4
Sanplatec PFA Test Tube Cap	<ul style="list-style-type: none"> ● Dedicated cap for Sanplatec test tube ● Made of PTFE ● It is knurled the side and good grip 		16249	For 15 and 20mL
Sanplatec PFA Test Tube Cap			16251	For 25mL
Sanplatec PFA Test Tube Cap			16252	For 50 and 60mL
Sanplatec PFA Test Tube	<ul style="list-style-type: none"> ● 2 line up: 20mL and 50mL ● It is translucent for easy view ● Comes with a dedicated PTFE cap 		16254	20mL-C (with Lid)
Sanplatec PFA Test Tube			16255	50mL-C (with Lid)
Acrylic Vacuum Desiccator (with drain valve)	<ul style="list-style-type: none"> ● Drain valves and adjustable vacuum gauge ● Draining is easy. ● vacuum gauge: 1pc (0~-0.1MPa) ● Intake and exhaust valve: 2pc ● Drain valve: 1pc, nozzle diameter: 1.5mm ● Gasket: silicon 		16532	UB-D type
Multi Wash Bottle	<ul style="list-style-type: none"> ● 360° washable ● The weight of glass ball makes it easy to wash ● Material: Main body: LDPE, Cap/Nozzle: PP, Tip nozzle: HDPE, weight: glass ● Individual package 		17000	500mL (Blue)
Multi Wash Bottle			17001	500mL (Red)
Multi Wash Bottle			17002	1L (Blue)
Multi Wash Bottle			17003	1L (Red)
EOG Sterile PP Bottle (wide mouth)	<ul style="list-style-type: none"> ● This bottle is sterilized by ethylene oxide ● Use as sample bottles for water analysis ● Individual package ● The cap stay removed when use 		17014c	100mL (200pcs)
EOG Sterile PP Bottle (wide mouth)			17015c	250mL (100pcs)
EOG Sterile PP Bottle (Wide Mouth)			17016c	500mL (100pcs)

EOG Sterile PP Bottle (Wide Mouth)			17017c	1L(50pcs)
EOG Sterile PP Bottle (wide mouth)	<ul style="list-style-type: none"> ● This bottle is sterilized by ethylene oxide gas (EOG). ● Use as sample bottles for water collection of water quality testing, food, cosmetics, chemicals, etc. ● With inner seal ● Individual package ● Both cap and inner seal stay removed when packed. 		17022c	100mL (200pcs)
EOG Sterile PP Bottle (wide mouth)			17023c	250mL (100pcs)
EOG Sterile PP Bottle (wide mouth)			17024c	500mL (100pcs)
EOG Sterile PP Bottle (wide mouth)			17025c	1L(50pcs)
PTFE Stir Bar (SA type)			17078	10mm(10pcs)
PTFE Stir Bar (SA type)			17079	15mm(10pcs)
PTFE Stir Bar (SA type)			17080	20mm(10pcs)
PTFE Stir Bar (SA type)			17081	25mm(10pcs)
PTFE Stir Bar (SA type)			17082	30mm(10pcs)
PTFE Stir Bar (SA type)			17083	35mm(10pcs)
PTFE Stir Bar (SA type)			17084	40mm(10pcs)
PTFE Stir Bar (SA type)			17085	45mm(10pcs)
PTFE Stir Bar (SA type)			17086	50mm(10pcs)
PTFE Stir Bar (SA type)			17087	55mm(10pcs)
Glove Box Compact	<ul style="list-style-type: none"> ● Compact size and light weight. ● Main body opening size : 250x325mm ● Valve: 2 pcs , Outer diameter of connection tube: ϕ 6mm, one-touch type 	● The outer dimension of main body includes the valves.	18003	GBC type
Sanplatec PFA Narrow Mouth Bottle (without inner seal)	<ul style="list-style-type: none"> ● Metal free. High quality PFA ● We conducted dissolution test ● Our factory conduct decompression test ● Leak-free with inner ring even after 1000 cycles ● Great airtightness ● Individual package 	● The total length includes a cap	18100	100mL
Sanplatec PFA Narrow Mouth Bottle (without inner seal)			18101	250mL
Sanplatec PFA Narrow Mouth Bottle (without inner seal)			18102	500mL
Sanplatec PFA Narrow Mouth Bottle (without inner seal)			18103	1000mL
Sanplatec PFA Narrow Mouth Bottle (with inner seal)	<ul style="list-style-type: none"> ● Metal free. High quality PFA ● We conducted dissolution test ● Our factory conduct decompression test ● Inner seal can be removed easily ● Great airtightness ● Individual package 	● The total length includes a cap	18104	100mL
Sanplatec PFA Narrow Mouth Bottle (with inner seal)			18105	250mL
Sanplatec PFA Narrow Mouth Bottle (with inner seal)			18106	500mL
Sanplatec PFA Narrow Mouth Bottle (with inner seal)			18107	1L
Sanplatec PFA Wide Mouth Bottle (with inner seal)	<ul style="list-style-type: none"> ● Metal free. High quality PFA ● We conducted dissolution test ● Our factory conduct decompression test ● Inner seal can be removed easily ● Great airtightness ● Individual package 	● The total length includes a cap	18112	100mL
Sanplatec PFA Wide Mouth Bottle (with inner seal)			18113	250mL
Sanplatec PFA Wide Mouth Bottle (with inner seal)			18114	500mL
Sanplatec PFA Wide Mouth Bottle (with inner seal)			18115	1000mL

Sanplatec Pre-Cleaned PFA Narrow Mouth Bottle (without inner seal)	●We clean Sanplatec PE Nar		18116	100mL(5pcs)
Sanplatec Pre-Cleaned PFA Narrow Mouth Bottle (without inner seal)			18117	250mL(5pcs)
Sanplatec Pre-Cleaned PFA Narrow Mouth Bottle (without inner seal)			18118	500mL(2pcs)
Sanplatec Pre-Cleaned PFA Narrow Mouth Bottle (without inner seal)			18119	1L
Sanplatec Pre-Cleaned PFA Wide Mouth Bottle (without inner seal)	●We clean Sanplatec PE Wide		18124	20mL(10pcs)
Sanplatec Pre-Cleaned PFA Wide Mouth Bottle (without inner seal)			18125	100mL(5pcs)
Sanplatec Pre-Cleaned PFA Wide Mouth Bottle (without inner seal)			18126	250mL(5pcs)
Sanplatec Pre-Cleaned PFA Wide Mouth Bottle (without inner seal)			18127	500mL(2pcs)
Sanplatec Pre-Cleaned PFA Wide Mouth Bottle (without inner seal)			18128	1L
Sanplatec PFA Narrow Mouth Wash Bottle	●This bottle can withstand an ●Main body,lid and nozzle are ●Autoclavable. ●Very useful in the field of mi	●The height does not include th	18133	100mL
Sanplatec PFA Narrow Mouth Wash Bottle			18134	250mL
Sanplatec PFA Narrow Mouth Wash Bottle			18135	500mL
Sanplatec PFA Narrow Mouth Wash Bottle			18136	1L
Sanplatec Pre-Cleaned PFA Narrow Mouth Wash Bottle	●We clean Sanplatec PE Narrow Mouth Bottle with pure water in clean rooms (cleanliness: class 100)and double packed in class 1000 environment.	●The height does not include the nozzle.	18137	100mL(5pcs)
Sanplatec Pre-Cleaned PFA Narrow Mouth Wash Bottle			18138	250mL(5pcs)
Sanplatec Pre-Cleaned PFA Narrow Mouth Wash Bottle			18139	500mL(2pcs)
Sanplatec Pre-Cleaned PFA Narrow Mouth Wash Bottle			18140	1L
Sanplatec PFA Wide Mouth Wash Bottle	●This bottle can withstand an ●Main body,lid and nozzle are ●Very useful in the field of mi	●The height does not include th	18141	100mL
Sanplatec PFA Wide Mouth Wash Bottle			18142	250mL
Sanplatec PFA Wide Mouth Wash Bottle			18143	500mL
Sanplatec PFA Wide Mouth Wash Bottle			18144	1L
Sanplatec Pre-Cleaned PFA Wide Mouth Wash Bottle	●We clean Sanplatec PE Wide	●The height does not include th	18145	100mL(5pcs)
Sanplatec Pre-Cleaned PFA Wide Mouth Wash Bottle			18146	250mL(5pcs)
Sanplatec Pre-Cleaned PFA Wide Mouth Wash Bottle			18147	500mL(2pcs)
Sanplatec Pre-Cleaned PFA Wide Mouth Wash Bottle			18148	1L
Sanplatec PFA Wide Mouth Gas Wash Bottle	●Excellent heat and chemical resistance. ●The outer diameter of tube x inner diameter(mm): 8x6 ●Filter ball pore size: 10 μ m ●Filter ball diameter: 25mm ●Filter ball materia: PTFE	●The height does not include the nozzle.	18149	100mL
Sanplatec PFA Wide Mouth Gas Wash Bottle			18150	250mL
Sanplatec PFA Wide Mouth Gas Wash Bottle			18151	500mL

Sanplatec PFA Wide Mouth Gas Wash Bottle			18152	1L
Sanplatec PFA Measuring Flask	<ul style="list-style-type: none"> ●Metal-free and high quality. ●Sanplatec original brand ●Chemical resistance (ex. str) ●The elution of impurities from ●It is produced by blow mould 		18153	25mL
Sanplatec PFA Measuring Flask			18154	50mL
Sanplatec PFA Measuring Flask			18155	100mL
Sanplatec PFA Measuring Flask			18156	250mL
Sanplatec PFA Measuring Flask			18157	500mL
PFA Tube(inch size)	<ul style="list-style-type: none"> ●Great non-adhesiveness and ●Suitable for electrical insulat ●Excellent weather tolerance. ●Translucent ●The maximum temperature is 		18490	1.59 × 3.17 × 0.79(10m)
PFA Tube(inch size)			18491	3.17 × 6.35 × 1.59(10m)
PFA Tube(inch size)			18492	4.35 × 6.35 × 1(10m)
PFA Tube(inch size)			18493	6.35 × 9.52 × 1.59(10m)
PFA Tube(inch size)			18494	7.52 × 9.52 × 1.59(10m)
PFA Tube(inch size)			18495	9.52 × 12.7 × 1(10m)
PFA Tube(inch size)			18496	10.7 × 12.7 × 1(10m)
PFA Tube(inch size)			18497	15.88 × 19.05 × 1.5(10m)
Sanplatec PTFE Beaker	<ul style="list-style-type: none"> ●According to our standards. ●High-quality beaker with exc ●Non-adhesiveness. Suitable 		19027	50mL
Sanplatec PTFE Beaker			19028	100mL
Sanplatec PTFE Beaker			19029	200mL
Sanplatec PTFE Beaker			19030	300mL
Sanplatec PTFE Beaker			19031	500mL
Sanplatec PTFE Beaker			19032	1000mL
Sanplatec PTFE Beaker			19033	2000mL
Large Size Desiccator SAWS type	<ul style="list-style-type: none"> ●Slide rail shelf type for putting ●Support: caster (with a stop) ●Flame/aluminium, caster / n ●Door opening size (mm): 497 ●Shelf size (mm): 460W × 460 ●Shelf material: SUS430 Pur ●Volume: 780L ●Shelf position: As you prefer ●Shelfuniform load: 12kg ●Silicagel container: 2pcs ●Total load: 120kg ●Main body weight: 63kg 	<ul style="list-style-type: none"> ●Valves can be attached via the (Upper and lower side: total 2 pc 	19604	SAWS-M-C
Large Size Desiccator SAWS type			19605	SAWS-M-CN
Large Size Desiccator SAWS type			19606	SAWS-M-CB
Large Size Desiccator SAWS type			19607	SAWS-H-CS
Large Size Desiccator SAWS type			19608	SAWS-H-CN
Large Size Desiccator SAWS type			19609	SAWS-H-CB
Large Size Desiccator SAWS type			19610	Shelf only (ST-1)
Large Size Desiccator SAWS type			19649	Shelf only (T-1)

Large Size Desiccator SAWD type (Auto dehumidification type)	<ul style="list-style-type: none"> ● Can be safely stored with a ● Humidity can be set (20% to ● Humidity control: Forced cor ● Dehumidification cycle: Repe ● Humidity setting: 20% to amb ● Humidity display: Humidity in ● Electric capacity: AC100V 1 ● Support: caster (with a stop ● Flame/aluminium, caster / ns ● Door opening size (mm):497 ● Shelf size (mm): 494W × 460 ● Shelf material: SUS430 Pur ● Volume:780L ● Shelf position: As you prefer ● Shelf uniform load: 12kg ● Total load: 120kg ● Main body weight: 63kg 	<ul style="list-style-type: none"> ● Valves can be attached via th (Upper and lower side: total 2 pc ● Shelf is not slide rail type. 	19623	SAWD-M-C
Large Size Desiccator SAWD type (Auto dehumidification type)			19624	SAWD-M-CN
Large Size Desiccator SAWD type (Auto dehumidification type)			19625	SAWD-M-CB
Large Size Desiccator SAWD type (Auto dehumidification type)			19626	SAWD-H-C
Large Size Desiccator SAWD type (Auto dehumidification type)			19627	SAWD-H-CN
Large Size Desiccator SAWD type (Auto dehumidification type)			19628	SAWD-H-CB
Large Size Desiccator SASD type(Auto dehumidification type)	<ul style="list-style-type: none"> ● Can be safely stored with a ● Humidity can be set (20% to ● Humidity control: Forced cor ● Dehumidification cycle: Repe ● Humidity setting: 20% to amb ● Humidity display: Humidity in ● Electric capacity: AC100V 1 ● Support: caster (with a stop ● Frame/aluminium, caster / ns ● Door opening size (mm): 497 ● Shelf size (mm): 494W × 460 ● Shelf material: SUS430 Pur ● Volume: 360L ● Shelf position: As you prefer ● Shelf uniform load: 12kg ● Total load: 120kg ● Main body weight: 36kg 	<ul style="list-style-type: none"> ● Valves can be attached via th (Upper and lower side: total 2 pc ● Shelf is not slide rail type. 	19629	SASD-M-C
Large Size Desiccator SASD type(Auto dehumidification type)			19630	SASD-M-CN
Large Size Desiccator SASD type(Auto dehumidification type)			19631	SASD-M-CB
Large Size Desiccator SASD type(Auto dehumidification type)			19632	SASD-P-C
Large Size Desiccator SASD type(Auto dehumidification type)			19633	SASD-P-CN
Large Size Desiccator SASD type(Auto dehumidification type)			19634	SASD-P-CB
Middle Size Desiccator SAG1 type	<ul style="list-style-type: none"> ● Material: main body/Acrylic ● Volume: 250L ● Door opening size (mm): 497 ● Shelf size (mm): 494W × 460 ● Shelf material: SUS430 Pur ● Volume: 360L ● Shelf position : As you prefer ● Shelf uniform load: 12kg ● Silicagel container: 1pc ● Total load: 120kg ● Main body weight: 36kg ● Caster with stoppers ● Total load: 120kg 	<ul style="list-style-type: none"> ● Valves can be attached via th (Upper and lower side: total 2 pc ● Shelf is not slide rail type. 	19655	SAG1-M-G
Middle Size Desiccator SAG1 type			19656	SAG1-M-GN
Middle Size Desiccator SAG1 type			19657	SAG1-M-GB
Middle Size Desiccator SAG1 type			19658	SAG1-M-C
Middle Size Desiccator SAG1 type			19659	SAG1-M-CN
Middle Size Desiccator SAG1 type			19660	SAG1-M-CB

Middle Size Desiccator SAG1 type			19661	SAG1-P-G
Middle Size Desiccator SAG1 type			19662	SAG1-P-GN
Middle Size Desiccator SAG1 type			19663	SAG1-P-GB
Middle Size Desiccator SAG1 type			19664	SAG1-P-C
Middle Size Desiccator SAG1 type			19665	SAG1-P-CN
Middle Size Desiccator SAG1 type			19666	SAG1-P-CB
Middle Size Desiccator SAG3 type	<ul style="list-style-type: none"> ●Material: main body/Acrylic ●Volume: 180L ●Door opening size (mm): 497 ●Shelf size (mm): 494W × 460 ●Shelf material: SUS430 Pur ●Volume: 360L ●Shelf position: As you prefer ●Shelf uniform load: 12kg ●Silicagel container: 1pc ●Total load: 120kg ●Main body weight: 20kg 	<ul style="list-style-type: none"> ●Valves can be attached via the (Upper and lower side; total 2 pc ●Shelf is not slide rail type. 	19667	SAG3-M-G
Middle Size Desiccator SAG3 type			19668	SAG3-M-GN
Middle Size Desiccator SAG3 type			19669	SAG3-M-GB
Middle Size Desiccator SAG3 type			19670	SAG3-P-G
Middle Size Desiccator SAG3 type			19671	SAG3-P-GN
Middle Size Desiccator SAG3 type			19672	SAG3-P-GB
Middle Size Desiccator SAG4 type		●Shelf is not slide rail type.	19650	Shelf only (T-2)
Middle Size Desiccator SAG4 type			19673	SAG4-M-G
Middle Size Desiccator SAG4 type			19674	SAG4-M-GN
Middle Size Desiccator SAG4 type			19675	SAG4-M-GB
Middle Size Desiccator SAG4 type			19676	SAG4-P-G
Middle Size Desiccator SAG4 type			19677	SAG4-P-GN
Middle Size Desiccator SAG4 type			19678	SAG4-P-GB
Middle Size Desiccator SAG5 type		●Shelf is not slide rail type.	19679	SAG5-M-G
Middle Size Desiccator SAG5 type			19680	SAG5-M-GN
Middle Size Desiccator SAG5 type			19681	SAG5-M-GB
Middle Size Desiccator SAG5 type			19682	SAG5-P-G
Middle Size Desiccator SAG5 type			19683	SAG5-P-GN
Middle Size Desiccator SAG5 type			19684	SAG5-P-GB
Middle Size Desiccator SAG1D type (auto dehumidification type)	<ul style="list-style-type: none"> ●This product can safely stor ●Humidity can be set (20% to ●Humidity control: Forced con ●Dehumidification cycle: Repe ●Humidity setting: 20% to amb ●Humidity display: Humidity in ●Electric capacity: AC100V 1 ●Material: main body/Acrylic ●Volume:250L ●Door opening size (mm): 497 ●Shelf size (mm): 494W × 460 ●Shelf material: SUS430 Pur ●Shelf position: As you prefer ●Shelf uniform load: 12kg ●Caster with stoppers ●Total load: 120kg ●Main body weight: 27kg 		19685	SAG1D-M-G
Middle Size Desiccator SAG1D type (auto dehumidification type)			19686	SAG1D-M-GN
Middle Size Desiccator SAG1D type (auto dehumidification type)			19687	SAG1D-M-GB
Middle Size Desiccator SAG1D type (auto dehumidification type)			19688	SAG1D-M-C
Middle Size Desiccator SAG1D type (auto dehumidification type)			19689	SAG1D-M-CN

Middle Size Desiccator SAG1D type (auto dehumidification type)			19690	SAG1D-M-CB
Middle Size Desiccator SAG1D type (auto dehumidification type)			19691	SAG1D-P-G
Middle Size Desiccator SAG1D type (auto dehumidification type)			19692	SAG1D-P-GN
Middle Size Desiccator SAG1D type (auto dehumidification type)			19693	SAG1D-P-GB
Middle Size Desiccator SAG1D type (auto dehumidification type)			19694	SAG1D-P-C
Middle Size Desiccator SAG1D type (auto dehumidification type)			19695	SAG1D-P-CN
Middle Size Desiccator SAG1D type (auto dehumidification type)			19696	SAG1D-P-CB
Portable Vacuum Sampling Box	<ul style="list-style-type: none"> ● This is an ultra-small size thin ● Valves are equipped as standard ● Check valves operates as reverse ● Light weight (approx.200g~600g) ● It is convenient not only for sampling ● Ball valve port: outer diameter: 10mm ● Check valve connection port: 1/8" ● Check valve operating pressure: 0.1MPa 	<ul style="list-style-type: none"> ● The total length includes the length of the sampling tube ● The outer dimension for round type 	19926	BOX Square type
Portable Vacuum Sampling Box			19927	BOX 50mL
Portable Vacuum Sampling Box			19928	BOX 100mL
Portable Vacuum Sampling Box			19929	BOX 150mL
Portable Vacuum Sampling Box			19930	BOX 300mL
Portable Vacuum Sampling Box			19931	BOX 500mL
Sanplatec Fluoro-Barrier PE Wide Mouth Botttle	<ul style="list-style-type: none"> ● The surface layer of the container is fluorine resin ● Fluorine is chemically bonded to the surface of the container (However, when the surface layer is damaged, the fluorine resin will be peeled off) ● The lot number of the fluorine resin is printed on the container ● Applications: for preserving cell culture, etc. ● With inner seal 	<ul style="list-style-type: none"> ● Since the substrate part is polycarbonate, it is not resistant to organic solvents ● Pay attention to handling the container ● Do not apply direct flame. ● Please clean before use. For cleaning, please refer to the manual. ● Each size is just for your reference 	20415c	100mL (200pcs)
Sanplatec Fluoro-Barrier PE Wide Mouth Botttle			20416c	250mL (100pcs)
Sanplatec Fluoro-Barrier PE Wide Mouth Botttle			20417c	500mL(50pcs)
Sanplatec Fluoro-Barrier PE Wide Mouth Botttle			20418c	1L(50pcs)
Sanplatec Fluoro-Barrier PE Wide Mouth Botttle			20419c	2L(30pcs)
PTFE Stirring Rod (with propeller)	<ul style="list-style-type: none"> ● Stirring rod with 5 blades ● The diameter of blades is small ● The shown size is for a blade diameter ● Material: PTFE (center core) ● The diameter of stainless core: 8mm ● The diameter of PTFE core: 10mm 		21163	25 φ
PTFE Stirring Rod (with propeller)			21164	35 φ
PTFE Stirring Rod (with propeller)			21165	50 φ
PTFE Stirring Rod (with propeller)	<ul style="list-style-type: none"> ● The blade propellers open when spinning ● The shown size is for a blade diameter ● Material: PTFE (center core) ● The diameter of stainless core: 8mm ● The diameter of PTFE core: 10mm 		21166	80 φ
PTFE Stirring Rod (with propeller)			21167	120 φ
PTFE Stirring Rod (with propeller)			21168	120 φ
PTFE Stirring Rod (with propeller)			21169	120 φ
PTFE Stirring Rod (with propeller)			21170	150 φ
PTFE Stirring Rod (with propeller)	<ul style="list-style-type: none"> ● The blade propellers open when spinning for high agitation . ● The shown size is for a blade diameter. (mm) ● The stainless core is also teflon-coating. ● Material :PTFE (center core: SUS) ● The diameter of PTFE core: 8mm 		21171	80 φ

PTFE Stirring Rod (with propeller)			21172	80 φ
PTFE Stirring Rod (with propeller)			21173	80 φ
PTFE Stirring Rod (with propeller)			21174	100 φ
PTFE Stirring Rod (with propeller)			21175	100 φ
PTFE Stirring Rod (with 3 propellers)	<ul style="list-style-type: none"> ● The 3-blade propellers enable ● The shown size is for a blade ● Material: PTFE (center core) ● The diameter of stainless core ● The diameter of PTFE core: 		21176	70 φ
PTFE Stirring Rod (with 3 propellers)			21177	70 φ
PTFE Stirring Rod (with 3 propellers)			21178	70 φ
PTFE Stirring Rod (with 2 propellers)	<ul style="list-style-type: none"> ● The 2-blade propellers enable ● The shown size is for a blade ● Material: PTFE (center core) ● The diameter of stainless core ● The diameter of PTFE core: 		21179	80 φ
PTFE Stirring Rod (with 2 propellers)			21180	80 φ
PTFE Stirring Rod (with 2 propellers)			21181	80 φ
Sanplatec PTFE Dish	<ul style="list-style-type: none"> ● Made of PTFE. Excellent chem ● Autoclavable ● Great non-adhesiveness ● A lid does not include. 		21221	60*50*15H
Sanplatec PTFE Dish			21222	90*80*15H
Sanplatec PTFE Dish			21223	100*90*20H
Sanplatec PTFE Dish			21224	120*110*20H
Sanplatec PTFE Dish			21225	150*140*20H
PFA Pipette	<ul style="list-style-type: none"> ● All made of PFA ● Useful application: Liquid tra ● Volume: 9mL 		21240	10 φ × 210L
Hydrofluoric Acid Washed PFA Bottle	<ul style="list-style-type: none"> ● Pack in a Class 100 clean ro ● It is the highest quality prod ● Double-packaged in unit 		23790	Narrow mouth 100mL without inner seal (5pcs)
Hydrofluoric Acid Washed PFA Bottle			23791	Narrow mouth 250mL without inner seal (5pcs)
Hydrofluoric Acid Washed PFA Bottle			23792	Wide mouth20mL without inner seal (10pcs)
Hydrofluoric Acid Washed PFA Bottle			23793	Wide mouth100mL without inner seal (5pcs)
Hydrofluoric Acid Washed PFA Bottle			23794	Wide mouth250mL without inner seal (5pcs)
Hydrochloric Acid Washed PFA Bottle	<ul style="list-style-type: none"> ● Pack in a Class 100 clean ro ● It is the highest quality prod ● Double-package in unit 		23795	Narrow mouth 100mL without inner seal (5pcs)
Hydrochloric Acid Washed PFA Bottle			23796	Narrow mouth 250mL without inner seal (5pcs)
Hydrochloric Acid Washed PFA Bottle			23797	Wide mouth 20mL without inner seal (10pcs)
Hydrochloric Acid Washed PFA Bottle			23798	Wide mouth 100mL without inner seal (5pcs)
Hydrochloric Acid Washed PFA Bottle			23799	wide mouth 250mL without inner seal (5pcs)
PFA Round Bottom Flask	<ul style="list-style-type: none"> ● This product is usually used a ● Do not use under decompress ● The main body of flask may be 		23898	100mL
PFA Round Bottom Flask			23899	250mL
PFA Round Bottom Flask			23900	500mL

Easy-Squeezey Wash Bottle	<ul style="list-style-type: none"> ●The soft material is used in ●Requires only minimal effort ●The bottle shape will recover ●Main body: LDPE, Cap, nozzle ●Individual package 		25000c	500mL(5pcs)
Gas Venting Wash Bottle	<ul style="list-style-type: none"> ●This item can eliminate the ●Adjust the inside pressure ●The nozzle tip is removal and ●Main body, Cap: LDPE, Noz 	●Make sure to tighten firmly th	25001	250mL
Gas Venting Wash Bottle			25002	500mL
Gas Venting Wash Bottle			25003	1L
Sanplatec Antistatic Wide Mouth Bottle	<ul style="list-style-type: none"> ●Semi-permanent antistat eff ●main body, Cap, Inner seal a ●The electrical resistance of ●It can reduce the loss of cor ●main body, Inner seal, Cap: F ●With inner seal 	<ul style="list-style-type: none"> ●The effect of this product is ●Do not apply to organic solve ●Each size is for reference onl 	25005c	100mL (200pcs)
Sanplatec Antistatic Wide Mouth Bottle			25006c	500mL (100pcs)
Sanplatec Antistatic Wide Mouth Bottle			25007c	1000mL (50pcs)
Sanplatec Antistatic Wide Mouth Bottle			25499c	250mL (100pcs)
Super Dry Keeper Extra (auto desiccator)	<ul style="list-style-type: none"> ●Tabletop type with digital re ●Humidity control is possible ●Timer operation of starting ●Dehumidifier: solid polymer e ●Dehumidification capacity: 4g ●The method of humidity conf ●The humidity control range: ●Humidity sensor: electrical c ●Humidity sensor accuracy: ± ●Power: single-phase AC100V ●Ambient service temperature ●Shelf can be adjusted to 15 ●Shelf: 3pcs 	●This product does not have a	25015	Extra
Super Dry Keeper Extra (auto desiccator)			00004	Shelf only(C-3)
Anesthesia Box (for small animals)	<ul style="list-style-type: none"> ●Anesthesia gas can be suppl ●Absorbent cottons can be pl ●You can also place a tray in ●Valve: 2pcs (nozzle diameter) ●Shelf: PVC punching board (25118	MAB-1
Anesthesia Box (for small animals)			25119	MAB-2
Glove Box Compact Vacuum Type	<ul style="list-style-type: none"> ●Excellent transparent view w ●The stand to put glove port ●Glove port I.D: 140mm ●Vacuum gauge range: 0~— ●Air-supply valve on main bo ●Air-exhaust valve on main b ●Air-exhaust valve on glove p ●Gasket material: silicon ●Glove material: natural rubbe ●piping tube: nylon, O.D 8mm 	<ul style="list-style-type: none"> ●Although this product is desig ●Vacuum gauge displays the di ●Do not operate inside during ●The outer dimension of the m 	25400	VG-C type
Sterilized Water Sampling Bottle with Hypo	<ul style="list-style-type: none"> ●This bottle includes sodium ●Suitable for bacteriological e ●A certain amount of sodium ●Main body: PE, Cap: PP ●Without Inner seal ●The cap is attached to the n ●Electronically Sterilized 		25425c	100mL Wide Mouth (200pcs)
Sterilized Water Sampling Bottle with Hypo			25426c	250mL Wide Mouth (100pcs)
Sterilized Water Sampling Bottle with Hypo			25427c	500mL Wide Mouth (100pcs)
Sterilized Water Sampling Bottle with Hypo			25428c	1L Wide Mouth (50pcs)
Sanplatec PFA Wide Mouth Bottle (without inner seal)	<ul style="list-style-type: none"> ●Metal free. High quality PFA ●We conducted dissolution te ●Our factory conduct decom ●Leak-free with inner ring eve ●Great airtightness ●Individual package 	●The total length includes a ca	13905	20mL
Sanplatec PFA Wide Mouth Bottle (without inner seal)			18108	100mL
Sanplatec PFA Wide Mouth Bottle (without inner seal)			18109	250mL
Sanplatec PFA Wide Mouth Bottle (without inner seal)			18110	500mL

Sanplatec PFA Wide Mouth Bottle (without inner seal)			18111	1000mL
Sanplatec PFA Wide Mouth Bottle (without inner seal)			25430	10mL
PTFE Stir Bar (SA type)			25464	10mm(50pcs)
PTFE Stir Bar (SA type)			25465	15mm(50pcs)
PTFE Stir Bar (SA type)			25466	20mm(50pcs)
PTFE Stir Bar (SA type)			25467	25mm(50pcs)
PTFE Stir Bar (SA type)			25468	30mm(50pcs)
PTFE Stir Bar (SA type)			25469	35mm(50pcs)
PTFE Stir Bar (SA type)			25470	40mm(50pcs)
Sanplatec Carboy with Nozzle (square bottle)			02502	20L
Sanplatec Carboy with Nozzle (square bottle)			25490	5L
Sanplatec Carboy with Nozzle (square bottle)			25491	10L
Sanplatec Fluoro-Barrier PE Tank	<ul style="list-style-type: none"> ● The surface layer of the container is fluorine resin. ● Fluorine is chemically bonded to the substrate. (However, when the surface layer is scratched, the fluorine resin will be removed.) ● The lot number of the fluorine resin is printed on the container. ● Applications: for preserving cell cultures, etc. ● With inner seal 	<ul style="list-style-type: none"> ● Since the substrate part is porous, please pay attention to handling the container. ● Do not apply direct flame. ● Please clean before use. ● For cleaning, use lukewarm water. ● Do not use cleaning brushes or detergents. ● Each size is for reference only. 	25500	20L Wide Mouth
Sanplatec Fluoro-Barrier PE Tank			25501	20L Narrow Mouth
Spray Bottle	<ul style="list-style-type: none"> ● 2 different mist types. ● No graduated scales 		25504	500mL
Drosophila Anesthesia Table	<ul style="list-style-type: none"> ● Carbon dioxide can be dyed. ● Anesthesia with CO₂ is less effective. ● Excellent operability with anesthesia. ● Bore size I.D: 8mm ● Top plate material: Plastic seal ● Anesthesia side dimension (φ): 100mm 	<ul style="list-style-type: none"> ● Do not clean or sanitize the anesthesia table. 	26010	KK-1
Butyl Long Gloves	<ul style="list-style-type: none"> ● The three-layered structure with a special synthetic rubber coating provides excellent resistance to aging in sunlight and ozone in the air. ● Excellent resistance to sunlight and atmospheric ozone. ● Three-layer construction: base/chloroprene(CR), middle/butyl(IIR), surface/chloroprene(CR) 		30220	7inch
Butyl Long Gloves			30221	8inch
UV Cut Container (Puratsubo)	<ul style="list-style-type: none"> ● 99.97% UV cut. Easy viewing ● The bottom is round to extra ● PP material. High chemical r ● Color: Main body: transparent 		26160	12mL(100pcs)
UV Cut Container (Puratsubo)			26161	24mL(100pcs)
UV Cut Container (Puratsubo)			26162	35mL(100pcs)
UV Cut Container (Puratsubo)			26163	58mL(100pcs)
UV Cut Container (Puratsubo)			26164	80mL(75pcs)
UV Cut Container (Puratsubo)			26165	125mL(50pcs)
Sterilized UV Cut Container (Puratsubo)	<ul style="list-style-type: none"> ● 99.97% UV cut. Easy viewing ● The bottom is round to extra ● PP material. High chemical r ● Color: Main body: transparent ● Sterilized with ethylene oxide 		26166	6mL(100pcs)
Sterilized UV Cut Container (Puratsubo)			26167	12mL(100pcs)
Sterilized UV Cut Container (Puratsubo)			26168	24mL(100pcs)
Sterilized UV Cut Container (Puratsubo)			26169	35mL(100pcs)
Sterilized UV Cut Container (Puratsubo)			26170	58mL(100pcs)
Sterilized UV Cut Container (Puratsubo)			26171	80mL(100pcs)
Sterilized UV Cut Container (Puratsubo)			26172	125mL(50pcs)

Sanplatec Fluoro-Barrier PE Wide Mouth Amber Bottle	<ul style="list-style-type: none"> ●The surface layer of the container is made of fluorine. ●Fluorine is chemically bonded to the surface. ●The lot number of the fluorine is recorded on the label. ●Applications: for preserving cells. ●With inner seal 	<ul style="list-style-type: none"> ●Since the substrate part is polycarbonate, please pay attention to handling the container. ●Do not apply direct flame. ●Please clean before use. ●For cleaning, use lukewarm water. Do not use cleaning brushes or detergents. ●Each size is for reference only. 	26230c	500mL(50pcs)
Sanplatec Fluoro-Barrier PE Wide Mouth Amber Bottle			26231c	1L(50pcs)
Sanplatec Fluoro-Barrier PE Wide Mouth Amber Bottle			26232c	2L(30pcs)
Ice Rack, Chill Pockets™	<ul style="list-style-type: none"> ●Main body part is made of stainless steel, which is easily heated and cool. The pockets to hold the tubing are made of aluminum. ●Tubes can fit well to the pockets to cool its whole body. ●Knurled type with whirl-stop function. You can open or close the tube cap with only one hand. ●The rack stabilizes on ice with moderate weight. ●Each aluminum pocket is independent. Ice can easily come into contact with each outer wall, making it easy to cool. ●The tube does not directly touch ice, reducing the risk of contamination ●Aluminum pockets fit 0.5ml, 1.5ml and 2.0ml of self-standing screw cap tube. ●When using an Eppen type tube, use a pocket non-knurled type(*There is no whirl-stop function) ●Knurled adaptable tubes Adapt to 0.5 ml, 1.5 ml, and 2.0 ml of assist tube, Watson tube, and SARSTEDT tube ●Material/main body: Stainless steel, Pocket: Aluminum, Joint Parts: PC ●Pocket placement pitch/23mm from center to center 		26559	6 Knurled Pockets
Ice Rack, Chill Pockets™			26560	12 Knurled Pockets
Ice Rack, Chill Pockets™			26561	24 Knurled Pockets
Ice Rack, Chill Pockets™			26562	6 Pockets
Ice Rack, Chill Pockets™			26563	12 Pockets
Ice Rack, Chill Pockets™			26564	24 Pockets
Gloveless Port	<ul style="list-style-type: none"> ●It can be attached to our glove box with a single touch. ●It can be used as a simple glove box when installed. ●The rubber material of the insertion part is NBR. 		26571	GRP-1(2pcs)
EB Sterilized Bottle	<ul style="list-style-type: none"> ●Sterilized by electron beam irradiation ●The main body and cap are individually packaged in a set. ●Special sealing structure provides high sealing performance even without inner seal. ●main body:PE、Cap:PP ●Without Inner seal 	<ul style="list-style-type: none"> ●The main body and cap may be slightly yellowed due to electron beam irradiation, but there is no problem in performance. ●Each size is just for your reference. ●The total length is with the cap attached. 	26583c	100mL (200pcs)
EB Sterilized Bottle			26584c	250mL (100pcs)
EB Sterilized Bottle			26585c	500mL (100pcs)
EB Sterilized Bottle			26586c	1L(50pcs)

Glove Box (Cubic type) Mobile lab space	<ul style="list-style-type: none"> ●The front opening is wide, making it easy to take things in and out. ●Easy to take on and off gloves. ●Lightweight and foldable for easy carrying. ●Foldable for storage. ●It can also be used as a simple gas purge. ●The supplied glove and replacement glove are simple disposable glove. ●Please use the Gas Purge glove (code 00193) for Gas Purge. 	<ul style="list-style-type: none"> ●The supplied glove and replacement glove are disposable glove. ●Please use the Gas Purge glove (code 00193) when performing Gas Purge. 	27001	MBSL-2
Screw Cap Opener	<ul style="list-style-type: none"> ●This screw cap opener allows you to open and close the cap without directly touching the tube. ●It allows you to easily open caps that are difficult to open after freezing. ●It can also be used in situations where you want to reduce the risk of contamination and avoid hazardous materials from sticking to your fingers. ●The tube can also be lifted for easy replacement in the rack. ●Press the knock button to release the tube. ●The tip cap is replaceable. ●Main body can be sprayed with alcohol ●Material/main body: PVC, back end: neodymium magnet, tip cap: polyurethane 	<ul style="list-style-type: none"> ●When using this product, please use our Ice lack (Code: 26559~26591). ●Not all screw cap tubes can be used with this product. 	26933	For 13~13.5mm (Cap)
Screw Cap Opener			26934	For 12~12.5mm (Cap)
Screw Cap Opener			26935	Tip cap only for 13 13.5mm (5pcs)
Screw Cap Opener			26936	Tip cap only for 12 12.5mm (5pcs)
Screw Cap Opener			27038	For 13.5 ~ 14mm
Screw Cap Opener			27039	Tip cap only for 13.5 14mm (5pcs)
Sanplatec Petri dish Multi-rack (Petri Multi-rack)	<ul style="list-style-type: none"> ●This is a petri dish for agar media. ●This rack can be used for all operations, such as loading, unloading, inverting, opening and closing the lid, carrying, and storing petri dishes. ●The top surface is deformable to widen the gap between the columns. ●It easy to carry stacked petri dishes. ●The rack can be inverted to hold petri dishes stably while they are in the rack. ●In the case of bacterial culture, the entire rack can be stored in an incubator. ●Autoclavable. ●Designed to minimize dead space when stored in an incubator or on a shelf. ●Material/main body: PP, washer: polycarbonate. 		26940	90 φ 5 Dish Model
Sanplatec Petri dish Multi-rack (Petri Multi-rack)			26941	90 φ 10 Dish Model

Formation Hood	<ul style="list-style-type: none"> ●It can be used as a table top hood by installing a fan unit on the ceiling, or as a table top clean bench by installing a clean unit. ●The ceiling has an opening for connection, so if you do not need the unit, please purchase a fitting for the formation hood (Code: 27902, 27903, 30112, 30113) and connect it. ●It can also be used for Table Top Clean Bench and Table Top Hood. ●To install the joint, simply fix it with screws. ●Material: Frame/Aluminum, Panel (except bottom)/PVC transparent board, Bottom/PVC white board 		26942	F-500H
Formation Hood			26943	F-700
Formation Hood			26944	F-1000
Formation Hood			26945	F-1000H
Fan Unit	<ul style="list-style-type: none"> ●Simply place it on top of the Formation Hood and tighten the included screws for easy installation. ●Power supply for 100V ●The power supply is a Japanese type outlet. ●When using the product overseas, please use a transformer or outlet adapter. 	<ul style="list-style-type: none"> ●The fan and metal parts are not resistant to corrosive chemicals. 	26947	1 Model
Fan Unit			26948	2 Model
Clean Unit	<ul style="list-style-type: none"> ●It can be used as a simple clean bench by simply placing it on top of the Formation Hood and fixing it with the included screws. ●Applicable main body: F-500H, F-700, F-1000, F-1000H ●Filter: HEPA filter, Non-woven fabric filter ●Airflow: 1m³/min ●Power: AC100V ●EOG gas sterilized. 		26953	1 Model
EOG Sterile Square Bottle	<ul style="list-style-type: none"> ●It is a square type, so there is no wasted space when storing. ●Can be used as water sampling bottles for water quality testing. ●All products are packaged individually. ●The inner seal and cap are both removed before packing. ●Material main body: PP, Cap: HDPE, Inner seal: LDPE ●With inner seal 	<ul style="list-style-type: none"> ●Each size is for reference only. 	26970	100mL (200pcs)
EOG Sterile Square Bottle			26971	250mL (100pcs)
EOG Sterile Square Bottle			26972	500mL (100pcs)
EOG Sterile Square Bottle			26973	1L(50pcs)
Sanplatec Round Container (Vessel)	<ul style="list-style-type: none"> ●The wide opening of the main body makes it ideal for storing and transporting paste materials and highly viscous samples. ●The lid has a handle to make it easy to carry. ●With inner seal ●Total length is the dimension when the lid is attached. ●Main body and inner sealMaterial: PE (polyethylene), CapMaterial: PP (polypropylene) 		26990	3.5L

Wash Bottle with Leak-proof venting cap (AFUREMA SENBIN)	<ul style="list-style-type: none"> ●The Auto Gas Vent is designed to prevent chemicals from dripping from the nozzle, even with highly volatile chemicals such as acetone and ethanol. ●When stored, the Auto Gas Vent lets the gas accumulated inside out through the ventilation holes on the side of the conical cap, preventing the internal pressure from rising excessively. ●When in use, the PP ball inside the conical cap rises due to a sudden increase in internal pressure caused by the gripping action of the hand. ●Main body, Cap: LDPE, Ball: PP, Pointy Cap: PE, Nozzle, Tip Nozzle, Ball Receptacle: HDPE 	<ul style="list-style-type: none"> ●If the bottle is tilted horizontally or downward, the liquid may not stop flowing or may seep out through the conical cap vents. ●If the product is tilted horizontally or downward, the liquid will not be able to stop and the chemical will seep out from the conical cap. ●When removing the conical cap to clean the inside, remove it slowly so as not to lose any small parts such as the ball or ball catcher. 	27023	250mL
Wash Bottle with Leak-proof venting cap (AFUREMA SENBIN)			27024	500mL
Fluoro-Barrier Wash Bottle (AFUREMA SENBIN)	<ul style="list-style-type: none"> ●The surface layer of the container (fluoro-barrier layer) is polyethylene fluorinated, which has excellent gas barrier (especially against organic solvents) and chemical resistance. ●Fluorine is chemically bonded, so the fluorinated surface layer has a permanent effect. (However, when the surface layer is scratched or shaved, polyethylene as the base material appears on the surface layer, and the effect of the fluorination treatment is eliminated. Please be careful when handling) ●The lot number of the fluorination treatment is shown on the product, and manufacturing control is performed. ●Main body, Cap: LDPE, Ball: PP, Conical Cap: PE, Nozzle, Tip Nozzle, Ball Cover: HDPE 	<ul style="list-style-type: none"> ●Since the substrate part is polyethylene, the heat resistance is about 60°C in the conventional way. ●Pay attention to handling the surface without scratches. ●Do not apply direct flame. ●Please clean before use. For cleaning, use lukewarm water or cleaning solution with little irritation. Do not use cleaning brushes or cleansers containing polishing agents that may cause scratches to the surface. ●Each size is just for your reference. 	27032	500mL
Wash Bottle (Shower type)	<ul style="list-style-type: none"> ●The Wash Bottle (Shower type) can greatly reduce the time and effort required for washing instruments. ●The diameter of the hole in the shower is highly accurate, ensuring stable injection. ●When stored, the plug cap prevents dust adhesion and contamination. ●Main body, Cap: LDPE, Over Cap: PP ●The product is packaged individually. 		27034	500mL
Wash Bottle (Shower type)			27035	1L
Adjustable Microplate Tilting Stand	<ul style="list-style-type: none"> ●It is easy to collect samples because the plate stand can be tilted freely. ●Pipette operation can be performed in a natural posture, reducing the burden. ●The plate stand can be tilted freely to make sample collection easier. 6-hole, 12-hole, 24-hole and 96-hole plates can be used. ●Effective dimensions of the plate stand (mm): 180 x 90 ●The maximum tilt angle is 35 degrees (the plate tilts only toward the front with a stopper, so there is no need to worry about it falling). ●Total length: 40mm 		27043	PA-1

Mouse Retention Glove (KAMARENAIINE)	<ul style="list-style-type: none"> ●Protects only the two fingers that are easily bitten by mouse. ●It is made of soft cowhide leather and is comfortable to wear. ●Soft cowhide leather for a comfortable fit. ●Material: Cowhide leather ●Cannot be washed in water or autoclaved. ●Can be sprayed with alcohol 		27281	SS-LH
Mouse Retention Glove (KAMARENAIINE)			27282	S-LH
Mouse Retention Glove (KAMARENAIINE)			27283	M-LH
Mouse Retention Glove (KAMARENAIINE)			27284	L-LH
Mouse Retention Glove (KAMARENAIINE)			27285	SS-RH
Mouse Retention Glove (KAMARENAIINE)			27286	S-RH
Mouse Retention Glove (KAMARENAIINE)			27287	M-RH
Mouse Retention Glove (KAMARENAIINE)			27288	L-RH
Rack with Tube Lock	<ul style="list-style-type: none"> ●Polycarbonate rack with anti-rotation function. ●Autoclavable. 		27298	12 Knurled Pockets
Rack with tube lock			27299	24 Knurled Pockets
Dish Holder for ϕ 35 mm (OKAMOCHI Series)	<ul style="list-style-type: none"> ●It is convenient for stacking and storing 35 ϕ dishes and moving multiple dishes. ●Space-saving design that does not take up much space in the incubator. ●The main body is made of stainless steel. ●The main body is made of stainless steel and can be autoclaved. ●Maximum number of stacks: 9, recommended 4-5. ●The higher the stacks, the less effective the upper dish is in preventing the liquid level from shaking. ●Suitable dish diameter (outer diameter of the lid): 39-40mm 		27300	35-DH
Medium Bottle Stand	<ul style="list-style-type: none"> ●This design prevents the pipette from touching the bottle opening. ●Two angle settings. ●Simple and compact design that can also be used inside a clean bench ●It can hold 500mL and 1l culture medium bottles at an angle, allowing pipette operation at an angle that is easy to work with. ●Alcohol spraying is possible. ●The square type medium bottle is designed so that the bottom corner is at the bottom, making pipetting easy even when the liquid volume is low. 		27304	BBS-1
Wash Bottle with Dimple	<ul style="list-style-type: none"> ●It has dimples to make it easier for people with small hands to grip. ●It has a bellows nozzle and a marble weight for multiple functions (can be used with various degrees of accuracy). ●It can be used from any angle. ●Main body: LDPE, Cap/nozzle: HDPE, Tube connector: PP, Tip nozzle: HDPE, Weight: GLASS, Weight cover: PP ●Products are packaged individually. 		27310	1L
Wash Bottle with Dimple			27310c	1L(5pcs)

PE Bottle with Dimple (KUBIRE Bottle)	<ul style="list-style-type: none"> ●The width of the dimple section is almost the same diameter as the width of a 500mL washing bottle, making it easy to hold. ●Universal design makes it easy to hold even 1L containers ●Main body:LDPE、Cap:HDPE 		27315	1L
PE Bottle with Dimple (KUBIRE Bottle)			27315c	1L(50pcs)
Inert Gas Circulation Purification Device	<ul style="list-style-type: none"> ●The system is connected to a glove box to circulate inert gas inside the box. ●The inert gas inside the box and moisture generated inside the box are removed by the adsorbent. ●The adsorbent columns are arranged in two cylinders in series. ●The adsorbent is replaceable. ●Circulating gas flow rate: 20L/min ●Power supply voltage: AC100V±10% 2A ●Moisture removal capacity: Dew point temperature -50°C td (*conditions apply) ●Moisture remover: Molecular Sieves 13X 		27360	SANWNO
Acrylic Vacuum Glove Box Clear View	<ul style="list-style-type: none"> ●Easy handling to move items in and out from the left side door. ●Use the right pass box to take in and out equipments during operation. ●Accessory Valve: body 2pcs, pass box 2pcs, glove port 1pc, pump connect side 2pcs, Vacuum gauge: body+pass box each 1pc (0~-0.1MPa) Glove+Glove band: each 1 pair Glove port door: 1pair (when decompressing) ●Valve nozzle diameter: φ9mm ●Internal dimension of pass box(mm):250Wx200Dx200H 	<ul style="list-style-type: none"> ●Do not operate with sudden decompression when using. Although this product is designed to withstand decompression in both of main body and pass box, it may cause failure. ●Vacuum gauge shows the pressure difference between ambient and inside of glove box. (Mpa) ●Do not operate when the inside is under decompression. ●The outer dimension of main body includes pass box and valves. 	27379	VG-600N
Stainless Vacuum Desiccator	<ul style="list-style-type: none"> ●Vacuum desiccator with excellent durability and airtightness ●The door is made of stainless steel frame and acrylic plate, so the contents can be checked. ●The SUS material eliminates the risk of corrosion when storing chemicals. ●Allowable vacuum: 0.1~101kPa ●Volume: Approx. 15L ●Vacuum gauge: 1 piece (for 0 to -0.1MPa) ●Vacuum gauge: 1 piece (for 0 to -0.1MPa) ●Valve: 2 pieces (nozzleouter diameterφ9 mm) ●Exterior: Stainless steel (SUS304) ●Interior: Stainless steel (SUS304) ●Shelf: Stainless steel, 1 piece (lower space H100mm, upper space H150mm) ●Door: Stainless steel frame (SUS304) + Acrylic plate 20tmm ●Door gasket: silicon rubber 		27400	SNV-250SUS

Butyl Long Gloves	<ul style="list-style-type: none"> ● Butyl structure has low gas permeability. (Suppresses gas transmission from workers' sweat, etc.) ● Suppresses outgas generation from glove ● For long time use, we recommend the flat hand type to reduce the workload. ● Total length(mm): 800 ● Material: Butyl rubber (IIR) 	27536	6inch 0.6t
Butyl Long Gloves		27537	8inch 0.4t
Butyl Long Gloves		27538	8inch 0.8t
Butyl Long Gloves		27539	7inch 0.4t
Small Container (Shading: type) (PURATSUBO)	<ul style="list-style-type: none"> ● It shuts out not only visible light, but also ultraviolet rays (200nm) to near-infrared rays (800nm). ● The bottom is round, so there is little residue. ● Made of PP (polypropylene), it has high chemical resistance. ● Main body color: Black ● Cap color: black 	27669	6mL(100pcs)
Small Container (Shading: type) (PURATSUBO)		27670	12mL(100pcs)
Small Container (Shading: type) (PURATSUBO)		27671	24mL(100pcs)
Small Container (Shading: type) (PURATSUBO)		27672	35mL(100pcs)
Small Container (Shading: type) (PURATSUBO)		27673	58mL(100pcs)
Small Container (Shading: type) (PURATSUBO)		27674	125mL(50pcs)
Small Container (Shading: type) (PURATSUBO)		27675	80mL(75pcs)
Mouse Retention Glove (Aramid Felt type) (KAMARENAIINE)	<ul style="list-style-type: none"> ● Protects the two fingers from being bitten by mouse. ● Protects your fingers from mouse teeth. ● For both left and right hand use ● It can be worn under the disposable glove to keep it clean. ● Material: Cowhide leather, aramid felt ● Cannot be washed in water or autoclaved. ● Can be sprayed with alcohol ● The fingers are reinforced with aramid felt. 	27875	SS
Mouse Retention Glove (Aramid Felt type) (KAMARENAIINE)		27876	S
Mouse Retention Glove (Aramid Felt type) (KAMARENAIINE)		27877	M
Mouse Retention Glove (Aramid Felt type) (KAMARENAIINE)		27878	L
Cryotube Dedicated Tongs (U type)	<ul style="list-style-type: none"> ● It is designed to follow the shape of cryotubes and Eppendorf tubes for secure gripping. ● The thin tip allows the user to wade through tubes in the tube rack and get a good grip. ● Made of stainless steel ● Autoclavable ● The tip is designed so that it does not touch the bottom of the desk or clean bench when placed. 	27879	U

Cryotube Dedicated Tongs (X type)	<ul style="list-style-type: none"> ● These tongs are specially designed for cryo tubes. ● The tongs are reverse-acting, so the tube can be held even when the hand is removed. ● It is designed to follow the shape of cryotubes and Eppendorf tubes for secure gripping. ● The thin tip allows the user to wade through tubes in the tube rack and get a good grip. ● Made of stainless steel ● Autoclavable ● The tip is designed so that it does not touch the bottom of the desk or clean bench when placed. ● With hook 		27880	X (Reverse effect type)
Cryotube Dedicated Tongs (L type)	<ul style="list-style-type: none"> ● These tongs have the same features as Tongs U (Code: 27879) and Tongs X (Code: 27880), and are specially designed for use with cryotubes and Eppendorf tubes. ● The tongs are designed for use with cryotubes and eppendorf tubes. ● Made of stainless steel ● Autoclavable ● The tip is designed not to touch the bottom of the desk or clean bench when placed. ● To prevent contamination, place the product so that the tip is facing upward. Be careful of the tip when doing so. 		28470	L type
Wash Bottle (Shower type)	<ul style="list-style-type: none"> ● The Wash Bottle (Shower type) can greatly reduce the time and effort required for washing instruments. ● The diameter of the hole in the shower is highly accurate, ensuring stable injection. ● When stored, the plug cap prevents dust adhesion and contamination. ● Main body, Cap: LDPE, Over Cap: PP ● The product is packaged individually. 		27881	1L
Well Plate Holder (OKAMOCHI Series)	<ul style="list-style-type: none"> ● It is convenient for stacking and storing well plate and moving multiple plates. ● The main body is made of stainless steel. ● The main body is made of stainless steel and can be autoclaved. ● Applicable plate size (mm): 86 x 128 ● * is the maximum height with the handle extended ● The main body is made of stainless steel and can be autoclaved. ● Approximate number of stacks: 2 to 3pcs(dish). ● The higher you stack them, the less effective the upper dish will be in preventing the liquid level from shaking. 		27883	WPH
Dish Holder (OKAMOCHI Series)	<ul style="list-style-type: none"> ● It is ideal for transferring liquid culture media that you do not want to shake. ● The main body is made of stainless steel and can be sterilized by autoclave. ● Approximate number of stacks: Maximum 9pcs, recommended 4-5pcs(dish) ● The higher you stack them, the less effective the upper dish will be in preventing the liquid level from shaking. ● * is the maximum height with the handle extended 		27885	60mm

Dish Holder (OKAMOCHI Series)		27886	90mm	
Dish Holder (OKAMOCHI Series)		27887	150mm	
Glove Box (Gas Purge type)	<ul style="list-style-type: none"> ● G and GS type are available. Main body and bottom part are detachable, allowing to put contents as much as possible. ● Bottom part of GS type is made of stainless. Worry-free when spilling chemicals inside. ● Pass box performs to decompress. ● The inner dimension of pass box (mm): 200W × 200D × 200H ● Stopcock for gas purge mounted on main body: 2 pcs (φ 6 tube of one touch type) ● Valve on pass box: 2 pcs (nozzle diameter φ 9mm) ● Vacuum gauge on pass box: 1 pc (for 0 ~ -0.1MPa) ● Glove and Glove band: 1 set ● Material of bottom board: acrylic (G type), stainless (GS type) 	<ul style="list-style-type: none"> ● Main body is not allowed to decompress. ● The outer dimension of main body includes pass box and valves. 	27889	G-SUS type
Glove Box (Aluminum frame type)	<ul style="list-style-type: none"> ● The main body is a lightweight glove box with an aluminum frame. ● All board parts are acrylic board. ※ Acrylic plate (white) at the bottom ● Two valves are provided for gas purge and dry air introduction. ● Valve main body: Needle valve x 2 (One-touch joint φ 6) ● Glove, glove stopper band: 1 pair of gloves, 2 O-rings 		27890	ALGB700 type
Glove Box (Aluminum frame type)			27891	ALGB1000 type
Glove Box (Aluminum frame type)			27892	ALGB1000H type
Nitric Acid Washed PE Bottle (Wide Mouth)	<ul style="list-style-type: none"> ● Nitric acid washed polyethylene bottles. ● Double-packed in pairs. ● Main body Material: PE, Cap Material: PP ● Without inner seal ● For more information, please contact us using the contact form. 		27954	100mL (5 pcs)
Nitric Acid Washed PE Bottle (Wide Mouth)			27955	500mL (1 pc)
Dish tray (φ 35mm)	<ul style="list-style-type: none"> ● It is easy to carry and stack 35mm dishes. ● Easy to carry and stack 35mm dishes in and out of incubator. ● It can be carried like a 6-well plate by simply placing a 35 mm dish on the tray. ● It is possible to exchange only some wells and to open and close the lid individually, which is not possible with the well plate. ● Alcohol disinfection is possible. 		27985	For 6 well
Beaker Lid (self-supporting type)	<ul style="list-style-type: none"> ● The lid is designed for use with TPX® beakers and PP beakers. ● The lid can stand on its own. ● A stopper on the back of the cap allows you to slide the lid to stabilize it in the slightly open position. ● Alcohol disinfection is possible. 		27986	1L
Beaker Lid (self-supporting type)			27987	2L
Beaker Lid (self-supporting type)			27988	3L
Beaker Lid (self-supporting type)			27989	5L

Acrylic Vacuum Desiccator (Package tester)	<ul style="list-style-type: none"> ● Large size rectangular package tester (vacuum desiccator). ● With liquid inside, you can suction without sucking the liquid into the pump. ● A drainage valve and an adjuster foot for height adjustment are included at the bottom for easy drainage. ● Also can be used as a vacuum desiccator. ● vacuum gauge: 1pc, 0~ -0.1 MPa ● Intake and exhaust valves: 2pcs, nozzle outer diameter 9mm ● Drain valve: 1pc, nozzle outer diameter 9mm ● Sample lift: 1 piece ● Packing: Silicon 		28235	OUB-D type
Glove Box (Gas Purge type)	<ul style="list-style-type: none"> ● This is a jumbo type of Gas Purge Glove Box DX-II (Code: 6590). ● The width of the main body part is 1200mm, which makes it possible to use a wide working space even when equipment is installed in the chamber. ● Accessories: glove 2pair, needle valve 4 pcs, vacuum gauge 1 pcs (pass box) ● Main body door opening dimensions: 1080 × 270(mm) ● It is designed to be easy to hold and can be moved delicately. ● The tip is flexible and can reach every corner of the container. ● Each bottle is individually wrapped for smooth opening. ● Sterilized by electron beam ● Size (mm)/Total length: 1500 ● Blade width: 13 ● Material/main body: Polycarbonate ● For more information about iP-TEC® series, please refer to the feature page on the website. 	<ul style="list-style-type: none"> ● The main body section cannot be depressurized. ● The main body outer dimension is the value including the pass box and valve. 	28238	DX-II JB
iP-TEC® Mini Scraper	<ul style="list-style-type: none"> ● It is designed to be easy to hold and can be moved delicately. ● The tip is flexible and can reach every corner of the container. ● Each bottle is individually wrapped for smooth opening. ● Sterilized by electron beam ● Size (mm)/Total length: 1500 ● Blade width: 13 ● Material/main body: Polycarbonate ● For more information about iP-TEC® series, please refer to the feature page on the website. 		28448	150-13 (1bag × 100)
Secondary Container of Live-Transport (for iP-TEC® flask-25)	<ul style="list-style-type: none"> ● Secondary container for iP-TEC® Flask-25 (Autoclavable) ● Two special trays are included for perfect-fitting. ● For more information about iP-TEC® series, please refer to the feature page on the website. ● 100% cellulose (code:28453) 		28451	1 Set
Secondary Container of Live-Transport (for iP-TEC® flask-25)			28452	Pallet of Live-transport
Secondary Container of Live-Transport (for iP-TEC® flask-25)			28453	Liquid Absorbing sheet

iP-TEC® Phase Change Material+36NS	<ul style="list-style-type: none"> ●Excellent temperature maintenance, keeps the temperature around 36°C for a long time. ●It is a heat storage material that can be used repeatedly. ●The container is tray-shaped and can be used as a heat retention tray for plates, etc. ●Alcohol spray OK ●Container Material: PVC ●On land: Designated combustible material under the Fire Service Law ●Marine: Non-hazardous material under the Ship Safety Law ●Air: Non-hazardous material under the Civil Aeronautics Law ●For more information about iP-TEC® series, please refer to the feature page on the website. 	<ul style="list-style-type: none"> ●Do not use autoclave or microwave oven. 	28457	1 pc
iP-TEC® Phase Change Material+24	<ul style="list-style-type: none"> ●It is a highly functional type of heat storage material that maintains a temperature of 24–25°C for a long time. ●It can be used when you want to transport the product in the same temperature range as the environmental temperature in a cell conditioning room (CPC) or medical site. ●It is a heat storage material that can be used repeatedly. ●The container is tray-shaped and can be used as a heat retention tray for plates, etc. ●Alcohol spray OK ●Container Material: PVC ●On land: Designated combustible material under the Fire Service Law ●Marine: Non-hazardous material under the Ship Safety Law ●Air: Non-hazardous material under the Civil Aeronautics Law ●For more information about iP-TEC® series, please refer to the feature page on the website. 	<ul style="list-style-type: none"> ●Do not use autoclave or microwave oven. 	28483	1 pc
iP-TEC® Temperature Stabilizer+5	<ul style="list-style-type: none"> ●High performance cold storage material that keeps samples near 5 to 6°C for a long time ●Stable stacking is possible. ●It can be used in the laboratory as a cold storage tray for plates, etc., in addition to being used for constant temperature transportation. ●Alcohol spray OK ●Material: PVC ●For more information about iP-TEC® series, please refer to the feature page on the website. 	<ul style="list-style-type: none"> ●Never use an autoclave or microwave oven. 	28648	1 pc
iP-TEC®PREMER BOX-V19	<ul style="list-style-type: none"> ●High performance insulated transport box using vacuum insulation material. ●For more information about iP-TEC® series, please refer to the feature page on the website. 	<ul style="list-style-type: none"> ●The secondary container and flask are sold separately. ●Performance will gradually degrade with age and storage environment. 	28460	BOX only

iP-TEC®STANDARD BOX-X13	<ul style="list-style-type: none"> ●Styrofoam insulation box with excellent cost performance ●Material: Interior: Aluminum evaporated foam PE, Packing: EPDM ●For more information about iP-TEC® series, please refer to the feature page on the website. 	<ul style="list-style-type: none"> ●The secondary container and flask are sold separately. 	28463	BOX only
iP-TEC® Light BOX-S6.6	<ul style="list-style-type: none"> ●The combination of the thermal insulation box and the phase change materials maintains the culture temperature range for a short time, making it ideal for short-time transportation between facilities. ●Weight: 1.5kg (with 2 phase change materials), 0.2kg (box only) ●For more information about iP-TEC® series, please refer to the feature page on the website. 	<ul style="list-style-type: none"> ●The secondary container and flask are sold separately. 	28466	Set (with 2 Temperature Stabilizer+36 NS)
iP-TEC® Light BOX-S6.6			28487	Set (with 2 Temperature Stabilizer+24)
iP-TEC® Tote BOX-6.6	<ul style="list-style-type: none"> ●This is a tertiary container that is easy to use for researchers to transport cells themselves. ●It looks like an ordinary bag, so it looks natural even when using public transportation. ●The weight is only about 3 kg even with the thermal storage material, making it convenient for transporting cells to research facilities within a short distance. ●Dimensions (mm): Box outer dimensions: 340 x 230 x H247, inner dimensions: 240 x 170 x H130 ●Material / Bag: Polyester ●Volume: 6.6l ●Weight / BOX only: approx. 830g 	<ul style="list-style-type: none"> ●This product does not include the heat storage material or secondary container. 	28543	BOX only
Rack with tube lock		<ul style="list-style-type: none"> ●Not autoclavable 	28469	Type B
Tube Freezing Holder	<ul style="list-style-type: none"> ●A maximum of six screw caps tubes can be set. ●The handle can be hooked on the edge of the liquid nitrogen container. ●The top lid is movable to prevent the tubes from floating up due to buoyancy. ●Some tubes may be slanted, in which case use the top cover to hold them in place. 	<ul style="list-style-type: none"> ●After freezing, the screws will freeze and become hard. Please wear cold-resistant gloves when opening the lid. ●The bottom of the tube must be self-supporting to be used with this product. 	28471	6 well
Adjustable Microplate Tilting Stand			28480	SUS
iP-TEC® PREMIER BOX-V8.5	<ul style="list-style-type: none"> ●High performance insulated transport box using vacuum insulation material. ●It is lightweight, compact, and easy to carry. ●With a shoulder belt for easy carrying. ●For more information about iP-TEC® series, please refer to the feature page on the website. 	<ul style="list-style-type: none"> ●The secondary container and flask are sold separately. ●The performance of the product will gradually decrease with age and storage environment. 	28482	BOX only

iP-TEC® Well Plate Cover	<ul style="list-style-type: none"> ● You can seal the plate without using the plate seal. ● The silicone cover is composed of an ultra-thin part that is necessary for ventilation and a thickness part that seals the container opening. ● The ultra-low hardness of the cover seals the container opening as if it were absorbing it. ● The cover is designed to reduce the amount of culture fluid during transportation and to prevent it from rippling, thus protecting the cells and tissues being transported. ● EOG sterilized ● Individual package ● Autoclavable ● Alcohol spray OK ● Size (mm)/140 × 97 × 10H (6, 12, 24well plate) /140 × 97 × 6H(96well plate) ● Material /Medical silicone ● For more information about iP-TEC® series, please refer to the feature page on the website. 	<ul style="list-style-type: none"> ● Not all dishes and plate brands are compatible. 	28488	for 6well(1 pc)
iP-TEC® Well Plate Cover			28489	for 6well(10 pcs)
iP-TEC® Well Plate Cover			28490	for 12well(1 pc)
iP-TEC® Well Plate Cover			28491	for 12well(10 pcs)
iP-TEC® Well Plate Cover			28492	for 24well(1 pc)
iP-TEC® Well Plate Cover			28493	for 24well(10 pcs)
iP-TEC® Well Plate Cover			28494	for 96well(1 pc)
iP-TEC® Well Plate Cover			28495	for 96well(10 pcs)
iP-TEC® Plank for Well Plate	<ul style="list-style-type: none"> ● Not sterilized ● Individual package ● Alcohol spray OK ● Size (mm)/128 × 85 × 5H(6, 12, 24, 96 well plate) ● Material/PVC ● For more information about iP-TEC® series, please refer to the feature page on the website. 	<ul style="list-style-type: none"> ● Not autoclavable ● Not all dishes and plate brands are compatible. 	28496	6well(1 pc)
iP-TEC® Plank for Well Plate			28497	6well(10 pcs)
iP-TEC® Plank for Well Plate			28498	12well(1 pc)
iP-TEC® Plank for Well Plate			28499	12well(10 pcs)
iP-TEC® Plank for Well Plate			28500	24well(1 pc)
iP-TEC® Plank for Well Plate			28501	24well(10 pcs)
iP-TEC® Plank for Well Plate			28502	96well(1 pc)
iP-TEC® Plank for Well Plate			28503	96well(10 pcs)

iP-TEC® Dish Cover	<ul style="list-style-type: none"> ● The silicone cover is composed of an ultra-thin part that is necessary for ventilation and a thickness part that seals the container opening. ● The ultra-low hardness of the cover seals the container opening as if it were absorbing it. ● The cover is designed to reduce the amount of culture fluid during transportation and to prevent it from rippling, thus protecting the cells and tissues being transported. ● EOG sterilized ● Individual package ● Autoclavable ● Alcohol spray OK ● Size (mm)/ ϕ 42.3 × 8H (ϕ 35mm)/ ϕ 62 × 9H (ϕ 60mm) / ϕ 95.5 × 11Hmm (ϕ 90mm) ● Material/ Medical silicone ● For more information about iP-TEC® series, please refer to the feature page on the website. 	<ul style="list-style-type: none"> ● Not all dishes and plate brands are compatible. 	28504	for ϕ 35mm(1 pc)
iP-TEC® Dish Cover			28505	for ϕ 35mm(1 pc)
iP-TEC® Dish Cover			28506	for ϕ 35mm(1 pc)
iP-TEC® Dish Cover			28507	for ϕ 35mm(1 pc)
iP-TEC® Dish Cover			28508	for ϕ 35mm(1 pc)
iP-TEC® Dish Cover			28509	for ϕ 35mm(1 pc)
iP-TEC® Plank for Dish Cover	<ul style="list-style-type: none"> ● Not sterilized ● Individual package ● Alcohol spray OK ● Size (mm)/ ϕ 43 × 5H (ϕ 35mm)/ ϕ 63 × 5Hmm (ϕ 60mm)/ ϕ 96 × 5H (ϕ 90mm) ● Material/PVC 	<ul style="list-style-type: none"> ● Not Autoclavable ● Not all dishes and plate brands are compatible. 	28510	ϕ 35mm(1 pc)
iP-TEC® Plank for Dish Cover			28511	ϕ 35mm(10 pcs)
iP-TEC® Plank for Dish Cover			28512	ϕ 60mm(1 pc)
iP-TEC® Plank for Dish Cover			28513	ϕ 60mm(10 pcs)
iP-TEC® Plank for Dish Cover			28514	ϕ 90mm(1 pc)
iP-TEC® Plank for Dish Cover			28515	ϕ 90mm(10 pcs)
iP-TEC® Secondary Container	<ul style="list-style-type: none"> ● The secondary container is made of polycarbonate for superior airtightness and strength. ● Shallow bottom design for easy loading and unloading of plates and dishes ● Designed to be stackable ● For more information about iP-TEC® series, please refer to the feature page on the website. 	<ul style="list-style-type: none"> ● When opening the secondary container, remove the metal fittings while holding the lid to prevent the contents from flying out. ● When autoclaving, separate the lid from the bottom of the container. 	28516	PC-0.5
Rack for Temperature Stabilizer	<ul style="list-style-type: none"> ● This is a convenient rack for iP-TEC® thermal storage material in an incubator. 		28522	For 8 pcs type
Rack for Temperature Stabilizer			28523	For 6 pcs type

iP-TEC® Mesh Cushion for Secondary Container	<ul style="list-style-type: none"> ●The reticular structure ensures high gas permeability and cushioning, as well as high holding capacity. ●Please wash and autoclave before use. ●Material: Elastomer ●20mm thickness: Height of the main body of the culture vessel: 10-19mm (not including the lid) ●30mm thickness: Height of the main body of the culture vessel: 18mm or more (height not including the lid) ●For more information about iP-TEC® series, please refer to the feature page on the website. 	28530	Thickness 20mm
iP-TEC® Mesh Cushion for Secondary Container		28531	Thickness 30mm
iP-TEC® Flask-25	<ul style="list-style-type: none"> ●This flask is designed to reduce the volume of culture medium when full, while the same culture area. ●The shape of the flask is designed for easy pipette operation and scraper operation. ●Stable stackability. ●Hydrophilic culture surface, sterilized by electron beam ●Plug Cap type ●Bags of 10 come with a convenient zipper. ●Volume / 30mL ●Effective cultivation area / 25cm² ●Material/Cap: PE, main body: polystyrene ●For more information about iP-TEC® series, please refer to the feature page on the website. 	28445	P-25 (10pcs × 10bags)
iP-TEC® Flask-25		28544	P-25 (10pcs × 1bag)
Sterilized Spoon (Electron beam sterilization)	<ul style="list-style-type: none"> ●The double packaging makes it ideal for use in clean rooms. ●Self-supporting design prevents spillage when placed on the table ●Slim head shape for easy insertion into the mouth of the container ●Pyrogen-free certification (0.05U/mL or lower) ●For powder and semi-solid sampling ●Material: Polypropylene ●Electron beam sterilization ●Packaging / Double individual packaging 	28611	1mL(100pcs)
Sterilized Spoon (Electron beam sterilization)		28612	5mL(100pcs)
Sterilized Spoon (Electron beam sterilization)		28613	10mL(100pcs)
Sterilized Spoon (Electron beam sterilization)		28614	20mL(100pcs)
Tweezers Stand	<ul style="list-style-type: none"> ●It can be used as a place to place instruments in a clean bench or on a laboratory table. ●Can be used in upward and downward directions according to the apparatus. 	28620	2 holes type
Tube Rack (Maintain the state of holding diagonally)	<ul style="list-style-type: none"> ●The transparent rack allows you to visually check the amount of chemicals inside. ●By sliding the middle plate, the tube or tube can be held at an angle. ●There is no need to angle the pipette by hand when operating it. ●Compact palm size ●Alcohol spraying is possible. 	28621	For 1.5mL tube

Tube Rack (Maintain the state of holding diagonally)			28622	15mL Centrifuge Tube
Tube Rack (Maintain the state of holding diagonally)			28623	50mL Centrifuge Tube
Barrier Pouch Barrier Pouch/ Water-absorbing sheet	<ul style="list-style-type: none"> ● It is a pressure-resistant, hermetically sealed pouch (secondary container) that can safely transport clinical and investigational specimens. ● No leakage even when internal pressure of 95kPa is applied in the temperature range of -40° C to +55° C. (Tested and proven by a third party public institution in Japan). ● Newly designed check valve system achieves high pressure resistance (patented). ● Barrier Pouch B4R holds one iP-TEC® secondary container (Code:28516). ● For more information about iP-TEC® series, please refer to the feature page on the website. 		28628	BRP-B4-R (20pcs)
Barrier Pouch Barrier Pouch/ Water-absorbing sheet			28630	BRP-ABS (20pcs) Water-absorbing sheet
Barrier Pouch Barrier Pouch/ Water-absorbing sheet			28651	BRB-A5-R (20pcs)
Beaker Rack	<ul style="list-style-type: none"> ● Coasters for preventing small beakers (20-50mL) from tipping over ● Supports beakers of different sizes at the top and bottom ● Can be used as an auxiliary device when pouring chemicals from heavy reagent bottles such as gallon bottles ● Code:28631 Main supported manufacturers: (beakers only): IWAKI 20mL, 50mL / HARIO 20mL, 50mL / SEKIYA 50mL ● Code: 28632 Main supported manufacturers: (beakers only): IWAKI 30mL, 50mL / HARIO 30mL, 50mL / SEKIYA 50mL 	<ul style="list-style-type: none"> ● Beaker sizes vary from manufacturer to manufacturer. ● Cannot be used with tall beakers or conical beakers. 	28631	For 20・50mL beaker
Beaker Rack			28632	For 30・50mL beaker
iP-TEC® Cell Transport Container #RPE12	<ul style="list-style-type: none"> ● This is a special container for transporting 12-well inserts. ● Closed (airtight) type and vented (gas exchange) type can be selected according to the application. ● Electron beam sterilized. ● Material/ Main body: PS, Cap: PE, Inner container: Medical silicone rubber ● For more information about iP-TEC® series, please refer to the feature page on the website. 	<ul style="list-style-type: none"> ● Culture inserts are not included in this product. 	28635	Closed type 12pcs (2/pcs × 6)
iP-TEC® Cell Transport Container #RPE12			28636	Vent type 12pcs (2/pcs × 6)

iP-TEC® Cell Transport Quadruple Container #24	<ul style="list-style-type: none"> ● It is transported by holding it firmly with the lid of the iP-TEC® secondary container and the mesh cushion. (The reticulation cushion is omitted in the photo.) ● The secondary container can hold up to 24 units (4 sets x 6 sets). ● The secondary container can hold up to 24 pieces (6 sets of 4 strands). ● Electronic ray sterilized (Code: 28637, 28638) ● Code 28639 is unsterilized. ● Alcohol spray OK (not autoclaveable). ● For more information about iP-TEC® series, please refer to the feature page on the website. 	<ul style="list-style-type: none"> ● Cap and main body are sold separately. ● The culture insert is not included in this product. 	28637	Silicon cap 60pcs (6pcs × 10)
iP-TEC® Cell Transport Quadruple Container #24			28638	Main Container 240pcs (24pcs/bag × 10)
iP-TEC® Cell Transport Quadruple Container #24			28639	Dedicated Rack (8 holes)
iP-TEC® Cell Sheet Transport Container	<ul style="list-style-type: none"> ● The cells are not shaken in the container by filling the container with culture medium and capping it without any air bubbles. ● The numbers in the table indicate the diameter of the storage sheet size. ● Electron beam sterilized ● Vent type ● Individually packaged (bag) ● Material/main body: polystyrene, cap: PE, medical silicone rubber ● Container size: 75 × 62mm, Total length: 20mm ● For more information about iP-TEC® series, please refer to the feature page on the website. 		28640	φ 38mm 6pcs (1pc/bag × 6)
iP-TEC® Cell Sheet Transport Container			28641	φ 50mm 6pcs (1pc/bag × 6)
iP-TEC® Flask-25 with vent cap	<ul style="list-style-type: none"> ● This flask is designed to reduce the volume of culture medium when full, while the same culture area. ● The shape of the flask is designed for easy pipette operation and scraper operation. ● Stable stackability. ● Hydrophilic culture surface, sterilized by electron beam ● Vent Cap type (with vent cap (made of silicon) to prevent air bubbles from entering) ● Bags of 5 come with a convenient zipper. ● Volume / 30mL ● Effective cultivation area / 25cm² ● Material/Cap: PE, main body: polystyrene ● For more information about iP-TEC® series, please refer to the feature page on the website. 		28642	50pcs (5pcs/bag × 10)
iP-TEC® Flask-25 with vent cap			28643	10pcs (10pcs/bag) (Vent cap only)

iP-TEC® PREMIER BOX-mini6.6	<ul style="list-style-type: none"> ●High performance insulated transport box using vacuum insulation material. ●It is lightweight, compact, and easy to carry. ●With a shoulder belt for easy carrying. ●A large transparent pocket on the top side fits an invoice (Japanese size). ●The logger can be visually checked from outside. Temperature and vibration loggers are neatly stored in the lid (recess size: 76 x 60 x Depth 33mm). ●For more information about iP-TEC® series, please refer to the feature page on the website. 	<ul style="list-style-type: none"> ●The secondary container and flask are sold separately. ●The temperature logger and vibration logger are sold separately. ●The performance will gradually deteriorate with age and storage environment. 	28645	BOX only
iP-TEC® Perfusion Culture Container with Ports	<ul style="list-style-type: none"> ●The inside and outside of the culture insert can be used for perfusion of culture medium. ●Electron beam sterilized. ●Material/main body: PS, Cap: PE, inner container: medical silicone rubber, inner seal with port: PC, compatible tubing: inner diameter 2.1-2.3 mm ●For more information about iP-TEC® series, please refer to the feature page on the website. 	<ul style="list-style-type: none"> ●Culture inserts are not included in this product. 	28653	Closed type (12pc)
iP-TEC® Perfusion Culture Container with Ports			28654	Vent type (12pcs)
Dropping Pipette (Square)			29196	2mL(500pcs)
Economy Dropping Pipette	<ul style="list-style-type: none"> ●Not sterilized 		29825	Non-sterilized (500pcs/case)
Economy Dropping Pipette			29826	Non-sterilized (500pcs/case)
Glove Box (Purge type)	<ul style="list-style-type: none"> ●GL type and GLS type glove boxes for Gas Purge are the larger models of G type (Code: 06903) and GS type (Code: 27889), and have the same specifications with the main body and removable bottom. ●The main body and bottom are removable. ●This model is equipped with a pass box that can reduce pressure as standard equipment, so it is possible to load and unload items in a more stable condition, allowing for speedy work. ●The GL type and GLS type are designed and manufactured for instruments. ●Internal dimensions of pass box (mm): 200W x 200D x 200H ●Faucet for main body Gas Purge: 2 pcs (φ6 tube one-touch type) ●Valve for pass box: 2pcs (nozzle diameter φ9mm) ●Vacuum gauge for pass box: 1 pc (for 0 to -0.1 MPa) ●Glove and glove band: 1 set 	<ul style="list-style-type: none"> ●The main body cannot be depressurized. 	30010	GL type
Glove Box (Purge type)			30011	GLS type

Acrylic Vacuum Glove Box Clear View M-750H type	<ul style="list-style-type: none"> ● This is a large type of M-600H (Code: 06593). It is possible to install a tall device in the glove box and work in an inert gas environment. ● Main body Acrylic (transparent): 1 vacuum gauge, 2 vacuum valves, 1 pair of gloves ● Pass box Acrylic (transparent): 1 vacuum gauge, 2 vacuum valves ● Vacuum valve: R1/4 nozzle 9mm 	<ul style="list-style-type: none"> ● Internal operation under reduced pressure is not possible. 	30012	M-750H
Glove Box (Aluminum frame type)	<ul style="list-style-type: none"> ● This glove box is designed to work with electronic balances in a glove box. ● The main body is a lightweight glove box with an aluminum frame. ● All the plates are made of acrylic. Only the bottom plate is made of stainless steel. ● The main body is equipped with two valves, which can be used for gas purge or dry air introduction. ● Valve main body: Needle valve x 2 (One-touch joint 6 φ) ● Pass box: Acrylic (transparent): 1 vacuum gauge, 2 vacuum valves ● Vacuum valve: R1/4 nozzle 9 φ mm ● Glove, glove stopper band: 1 pair of glove, 2 O-rings 		30013	WG-ALGB
O-Ring for Glove Port			30230	O-ring for 7 inch
O-Ring for Glove Port			30231	O-ring for 8 inch
Formation Hood (Assembly type)	<ul style="list-style-type: none"> ● Assemblable without tools. ● It can be used as a table top hood when the fan unit is installed on the top surface, or as a table top clean bench when the clean unit is installed. ● When not in use, it can be disassembled for storage (special trunk for storage is sold separately). ● Safe and convenient: The front door is equipped with a constant-load spring, which allows you to stop the front door at any position. 	<ul style="list-style-type: none"> ● This product does not have a bottom plate. ● The special storage trunk, fan unit, and clean unit are not included. (They are sold separately.) 	30301	BF-500H
Formation Hood (Assembly type)			30302	BF-700
Formation Hood (Assembly type)			30303	BF-1000
Formation Hood (Assembly type)			30304	BF-1000H
Formation Hood (Assembly type)			30305	BF-500H (Dedicated Trunk)
Formation Hood (Assembly type)			30306	BF-700 (Dedicated Trunk)
Formation Hood (Assembly type)			30307	BF-1000 (Dedicated Trunk)
Formation Hood (Assembly type)			30308	BF-1000H (Dedicated Trunk)

ATTRAN®	<ul style="list-style-type: none"> ● This product has been created by applying special polymer on the polylactic acid non-woven fabric. ● Thus it is uniquely flexible and high in tracking to a cultured cell sheet and the affected body part. ● Safety 1) Primary skin irritation tested 2) Cytotoxicity evaluated 3) Gamma-ray sterilized (sterility test confirmed) 4) Bacterial endotoxins test performed ● Sterilized by gamma radiation 	<ul style="list-style-type: none"> ● This product is a laboratory device, not a medical device or medical instrument. ● Do not use on humans for any purpose, such as manufacturing, quality control, diagnosis, treatment, and research. 	30501	ATNC-030(10 pieces/bag × 5)
ATTRAN®			30502	ATNS-7686 (Pack of 10 pieces)
Sanplatec PFA Bottle GL45 type	<ul style="list-style-type: none"> ● The GL45 port is compatible with other manufacturers' GL45 caps. ● Cap with GL45 port can be used for loading and unloading of chemicals in a closed system. ● The main body is thickness and strong. ● The 1L type has two dimples on the side of the main body to make it easier to hold. ● Heat and cold resistance (-196° C to 260° C) ● Chemical resistance to acids, alkalis and organic solvents ● Metal-free, weather-resistant, low adsorption. ● Without inner seal 	<ul style="list-style-type: none"> ● Total length indicates the dimension when Cap is installed. 	30703	500mL
Sanplatec PFA Bottle GL45 type			30704	1L
Sanplatec Pyrogen-free PFA Narrow Mouth Bottle without inner seal	<ul style="list-style-type: none"> ● This is a new cleaning service product that removes pyrogen from objects. ● These PFA bottles have been cleaned with pyrogen-free water (PF water) in a clean room environment. ● The endotoxin competence level of PF water is 0.002 EU/mL. ● The PF water removes not only pyrogen but also particles from the object. ● After PF washing, the bottles are double-packed (in units of the number of bottles) in a high-level clean environment. ● Without inner seal 		30705	100mL(5pcs)
Sanplatec Pyrogen-free PFA Narrow Mouth Bottle without inner seal			30706	250mL(5pcs)
Sanplatec Pyrogen-free PFA Narrow Mouth Bottle without inner seal			30707	500mL(2pcs)
Sanplatec Pyrogen-free PFA Narrow Mouth Bottle without inner seal			30708	1L(1pc)

Sanplatec Pyrogen-free PFA Wide Mouth Bottle without inner seal	<ul style="list-style-type: none"> ● This is a new cleaning service product that removes pyrogen from objects. ● These PFA bottles have been cleaned with pyrogen-free water (PF water) in a clean room environment. ● The endotoxin competence level of PF water is 0.002 EU/mL. ● The PF water removes not only pyrogen but also particles from the object. ● After PF washing, the bottles are double-packed (in units of the number of bottles) in a high-level clean environment. ● Without inner seal 		30709	100mL(5pcs)
Sanplatec Pyrogen-free PFA Wide Mouth Bottle without inner seal			30710	250mL(5pcs)
Sanplatec Pyrogen-free PFA Wide Mouth Bottle without inner seal			30711	500mL(2pcs)
Sanplatec Pyrogen-free PFA Wide Mouth Bottle without inner seal			30712	1L(1pc)
Sanplatec Bio Plastic PE Wide Mouth Bottle	<ul style="list-style-type: none"> ● This bottle is made by blending ordinary polyethylene with bio-polyethylene, which reduces CO2 emissions by up to 70%(※1). ● The performance is the same as conventional petroleum-derived polyethylene. ● Unlike biodegradable plastic, it does not decompose naturally. ● It can be incinerated as waste with zero CO2 emissions. ● Bio-polyethylene is used for the main body, inner seal and cap. (※1)The main raw material, sugarcane, absorbs CO2 through photosynthesis during its growth stage, so when it is incinerated as waste, CO2 emissions can be regarded as zero (carbon neutral). Therefore, even including the manufacturing and transportation processes, CO2 emissions can be reduced by up to 70% compared to conventional petroleum-derived polyethylene, which is effective from the perspective of corporate social responsibility by breaking away from dependence on petroleum and contributing to the prevention of global warming. 	<ul style="list-style-type: none"> ● Each size is just for your reference. ● The total length is with the cap attached. 	30804c	100mL (200pcs)
Sanplatec Bio Plastic PE Wide Mouth Bottle			30805c	250mL (100pcs)
Sanplatec Bio Plastic PE Wide Mouth Bottle			30806c	500mL (100pcs)
Sanplatec Bio Plastic PE Wide Mouth Bottle			30807c	1L(50pcs)
Sanplatec PE Tough Bottle	<ul style="list-style-type: none"> ● This is a thicker type of liethylene bottle compared to our conventional products. ● Material:main body、Cap/HDPE、Inner seal/LDPE ● With inner seal 	<ul style="list-style-type: none"> ● Each size is just for your reference. ● Total length is with the cap attached. 	30816c	100mL (200pcs)
Sanplatec PE Tough Bottle			30817c	250mL (100pcs)
Sanplatec PE Tough Bottle			30818c	500mL (100pcs)
Sanplatec PE Tough Bottle			30819c	1L(50pcs)

Sanplatec Antistatic Wash Bottle	<ul style="list-style-type: none"> ●The main body, cap and nozzle all contain antistatic agents. ●The surface resistance (Ω) of the main body is 108-1011 on average. ●The special polymer-type antistatic agent provides a semi-permanent antistatic effect. ●Material: Cap, nozzle/HDPE, main body/LDPE 	<ul style="list-style-type: none"> ●The effect of this product is antistatic. There is no static elimination effect of an object that is already charged. ●Do not apply to organic solvents. In addition, acid and alkali may degrade the performance. ●Each size is for reference only. 	30826	500mL
Sanplatec Antistatic Wash Bottle			30826c	500mL (100pcs)
Sanplatec Pyrogen-free PE Wide Mouth Bottle	<ul style="list-style-type: none"> ●This is a new cleaning service product that removes pyrogen from objects. ●These PE bottles have been cleaned with pyrogen-free water (PF water) in a clean room environment. ●The endotoxin competence level of PF water is 0.002 EU/mL. ●The PF water removes not only pyrogen but also particles from the object. ●After PF washing, the bottles are double-packed (in units of the number of bottles) in a high-level clean environment. ●With inner seal 		30828	100mL(5pcs)
Sanplatec Pyrogen-free PE Wide Mouth Bottle			30829	250mL(5pcs)
Sanplatec Pyrogen-free PE Wide Mouth Bottle			30830	500mL(2pcs)
Sanplatec Pyrogen-free PE Wide Mouth Bottle			30831	1L(1pc)
Sanplatec Pyrogen-free PE Wide Mouth Bottle			30832	2L(1pc)
ONE FINGER Spray Bottle	<ul style="list-style-type: none"> ●Compared to conventional lever-type spray bottles, this product can be dispensed with just one thumb, making it easy to use continuously. 		30905	500mL (Black)
ONE FINGER Spray Bottle			30906	500mL (White)
Sanplatec Bio Plastic PE Disposable Cup (PE Beaker)	<ul style="list-style-type: none"> ●This Cup is made by blending ordinary polyethylene with bio-polyethylene, which reduces CO2 emissions by up to 70%(※1). ●The performance is the same as conventional petroleum-derived polyethylene. ●Unlike biodegradable plastic, it does not decompose naturally. ●It can be incinerated as waste with zero CO2 emissions. (※1)The main raw material, sugarcane, absorbs CO2 through photosynthesis during its growth stage, so when it is incinerated as waste, CO2 emissions can be regarded as zero (carbon neutral). Therefore, even including the manufacturing and transportation processes, CO2 emissions can be reduced by up to 70% compared to conventional petroleum-derived polyethylene, which is effective from the perspective of corporate social responsibility by breaking away from dependence on petroleum and contributing to the prevention of global warming. 	<ul style="list-style-type: none"> ●Each size is just for your reference. 	31030c	100mL (1000pcs)
Sanplatec Bio Plastic PE Disposable Cup (PE Beaker)			31032c	500mL (250pcs)
Sanplatec Bio Plastic PE Disposable Cup (PE Beaker)			31033c	2L(100pcs)

Sterilized Scoop (Electron beam sterilization)	<ul style="list-style-type: none"> ● A disposable scoop that can be placed directly on the balance after weighing. ● Self-supporting shape to prevent spillage when placed ● Slim head shape for easy insertion into the mouth of the container ● Pyrogen-free certification (0.05U/mL or lower) ● Suitable for powder and semi-solid sampling. ● Material: Polypropylene ● Sterilized by electron beam ● Packaging / Double individual packaging 		31040	100mL(50pcs)
Sanplatec Antistatic PE Disposable Cup (PE Beaker)	<ul style="list-style-type: none"> ● Antistatic function is added to the general PE Disposable Cup (PE Beaker). ● It is effective in preventing the adhesion of particles and the risk of insulation breakdown of contents. ● The special polymer-type antistatic agent is expected to have a semi-permanent antistatic effect. ● This product eliminates the hassle of sticking to walls when handling powders. ● Main body: Surface resistivity (Ω): 10^8 to 10^{11} Refer to JIS K 6911 General test method for thermosetting plastics (5.13 resistivity) 	<ul style="list-style-type: none"> ● The effect of this product is antistatic. There is no static elimination effect of an object that is already charged. ● Do not apply to organic solvents. In addition, acid and alkali may degrade the performance. ● Each size is for reference only. 	31041c	100mL (1000pcs)
Sanplatec Antistatic PE Disposable Cup (PE Beaker)			31042c	500mL (250pcs)
Sanplatec Antistatic PE Disposable Cup (PE Beaker)			31043c	1L (100pcs)
Tube Rack (For medium bottle only)	<ul style="list-style-type: none"> ● It is a tube holder that attaches to the mouth of a 500mL reagent bottle. ● Most suitable for use with 50mL centrifuge tubes, but can also be used with 15mL types ● By attaching directly to the bottle, mistakes such as misplacement can be prevented and work can be carried out. ● Directly insert the tube into the bottle and dispense the solution. 		31201	S
Spoon (Made of silicon)	<ul style="list-style-type: none"> ● It is very useful for removing highly viscous contents from containers. ● It is effective in reducing liquid residue. ● Heat resistant up to 200°C. ● Autoclavable. 		59000	S
Spoon (Made of silicon)			59001	M
Spoon (Made of silicon)			59002	L
Cushion for Protecting Glass Bottles (safety netron)	<ul style="list-style-type: none"> ● Perfect for protection of reagent bottles, industrial chemicals bottles, and various glass devices ● Stretch material can fit any size of containers. ● Excellent cold resistance, no change in elasticity when placed in a freezer. ● Easy cutting with scissors and cutters ● Material: Ethylene-vinyl acetate copolymer 		80900	DN-1 100M
Cushion for Protecting Glass Bottles (safety netron)			80901	DN-2 100M
Cushion for Protecting Glass Bottles (safety netron)			80902	DN-3 50M
Cushion for Protecting Glass Bottles (safety netron)			80903	DN-4 50M

Cushion for Protecting Glass Bottles (safety netron)			80904	DN-5 50M
Cushion for Protecting Glass Bottles (safety netron)			80905	DN-6 50M
Color Nozzle for Wash Bottle Replacement	<Compatible Wash bottles> ●Sanplatec® Round Wash Bottle (Code: 02117-02119), Square Wash Bottle (Code: 02115, 02116), Multi Wash Bottle (Code: 17000-17003), Wash Bottle (Softy) (Code: 25000), Wash Bottle (27023, 27024), Fluorobarrier Wash Bottle (27032) ●Colors: red, green, yellow, white, and blue.		90826	Red
Color Nozzle for Wash Bottle Replacement			90827	Green
Color Nozzle for Wash Bottle Replacement			90828	Yellow
Color Nozzle for Wash Bottle Replacement			90829	White
Color Nozzle for Wash Bottle Replacement			90830	Blue
Sanplatec Carboy with Two Ports	●Two-ports polyethylene carboy.		91530	10L
Sanplatec Carboy with Two Ports			91531	20L